



## Invasive Blood pressure Cables and transducers




### Instrument Connector



### Yoke type



### transducers connector



### Type Disignation

X □ □ □ □ □

Yoke type: A: Appott 6pin; B: Utah 4pin; C: Edward 5pin; D: BD 7pin; E: Medix 6pin; F: Philips 4pin

Compatible Brand code, as 0014 indicated Marquette, 0015 indicated Spacelabs.

IX: Indicated IBP Cable and Pressure Transducers

## Invasive Blood pressure Cables and transducers

| Compatible Brand | Compatible Models No.                                                                     | OEM P/N | MED-LINK REF NO. | Description (Name,Length,Plug)                                         | INSTRUMENT CONNECTOR                                                                      | Yoke (transducers side)                                                               |
|------------------|-------------------------------------------------------------------------------------------|---------|------------------|------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Chenwei</b>   |                                                                                           |         | EX0017A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>6p6c jack, |  251   |    |
|                  |                                                                                           |         | EX0017B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>Din 2.0 4pin |  251   |    |
|                  |                                                                                           |         | EX0017C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>5pin,      |  251   |    |
|                  |                                                                                           |         | EX0017D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>7p jack,       |  251   |    |
|                  |                                                                                           |         | EX0017E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>5pin        |  251   |    |
|                  |                                                                                           |         | EX0017F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Readet 4pin>Round 4p  |  251   |    |
| <b>Datex</b>     |  X0030  |         | EX0030A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>6p6c jack,   | 252    |    |
|                  |                                                                                           |         | EX0030B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>Din 2.0 4pin   | 252  |  |
|                  |                                                                                           |         | EX0030C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>5pin,        | 252  |  |
|                  |                                                                                           |         | EX0030D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>7p jack,         | 252  |  |
|                  |                                                                                           |         | EX0030E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>5pin          | 252  |  |
|                  |                                                                                           |         | EX0030F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>Round 4p    | 252  |  |
| <b>Fukuda</b>    |  X0026 |         | EX0030G          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack, | 253  |  |
|                  |                                                                                           |         | EX0026B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin | 253  |  |
|                  |                                                                                           |         | EX0026           | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin,      | 253  |  |
|                  |                                                                                           |         | EX0026D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,       | 253  |  |
|                  |                                                                                           |         | EX0026E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin        | 253  |  |
|                  |                                                                                           |         | EX0026F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p  | 253  |  |



## Invasive Blood pressure Cables and transducers

| Compatible Brand    | Compatible Models No.                                                                                                | OEM P/N    | MED-LINK REF NO. | Description (Name,Length,Plug)                                             | INSTRUMENT CONNECTOR                                                                      | Yoke ( transducers side)                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------|------------|------------------|----------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>GE-marquette</b> | GE Solar, Dash, Eagle, Tram<br><br> |            | EX0014A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>6p6c jack, | 413    |    |
|                     |                                                                                                                      | 700078-001 | EX0014B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11p>Din2.0 4p      | 413    |    |
|                     |                                                                                                                      |            | EX0014C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>5pin,      | 413    |    |
|                     |                                                                                                                      |            | EX0014D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>7p jack,       | 413    |    |
|                     |                                                                                                                      |            | EX0014E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>5pin        | 413    |    |
|                     |                                                                                                                      |            | EX0014F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>Round 4p  | 413    |    |
| <b>Lohmeier</b>     | M010<br><br>                      |            | EX0019A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>6p6c jack,   | 079    |    |
|                     |                                                                                                                      |            | EX0019B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>Din 2.0 4pin   | 079  |  |
|                     |                                                                                                                      |            | EX0019C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>5pin,        | 079  |  |
|                     |                                                                                                                      |            | EX0019D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>7p jack,         | 079  |  |
|                     |                                                                                                                      |            | EX0019E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>5pin          | 079  |  |
|                     |                                                                                                                      |            | EX0019F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>Round 4p    | 079  |  |
| <b>Mindray</b>      | PM6000,6800, T5, T8<br><br>       |            | EX0016A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack,     | 519  |  |
|                     |                                                                                                                      |            | EX0016B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin     | 519  |  |
|                     |                                                                                                                      |            | EX0016C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin,          | 519  |  |
|                     |                                                                                                                      |            | EX0016D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,           | 519  |  |
|                     |                                                                                                                      |            | EX0016E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin            | 519  |  |
|                     |                                                                                                                      |            | EX0016F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p      | 519  |  |

## Invasive Blood pressure Cables and transducers

| Compatible Brand                                                                    | Compatible Models No.                                                                                                                                                                                                                                                                                                                                 | OEM P/N     | MED-LINK REF NO. | Description (Name,Length,Plug)                                         | INSTRUMENT CONNECTOR                                                                      | Yoke ( transducers side)                                                              |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|------------------|------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Spacelabs<br>Mindray<br>Goldway                                                     | Reusable and non-sterile. It may be gas sterilized. (Datex-Ohmeda TC-VTK.) For use with 690-0008-00, 690-0009-00, 690-0010-00, 690-0011-00, 690-0014-01, and 690-0015-01 transducers and all Spacelabs Medical monitoring systems. Use on 90496 also requires adapter cable 700-0028-00, sold separately.<br><u>Mindray</u> : PM7000, 8000,9000,6201. |             | EX0015A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6p6c jack,    |  377   |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0015B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Din 2.0 4pin    |  377   |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0015C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>5pin,         |  377   |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       | 690-0021-00 | EX0015D          | BD Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>7p jack,          |  377   |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0015E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6pin           |  377   |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0015F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Round 4p     |  377   |    |
|    |                                                                                                                                                                                                                                                                                                                                                       |             |                  |                                                                        |                                                                                           |                                                                                       |
| Philips                                                                             | M1567A,M1568A kits                                                                                                                                                                                                                                                                                                                                    |             | EX0018A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack, | 255    |    |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0018B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin | 255  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0018C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin,      | 255  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0018D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,       | 255  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0018E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>5pin        | 255  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       | M1634A      | EX0018F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p  | 255  |  |
|  |                                                                                                                                                                                                                                                                                                                                                       |             |                  |                                                                        |                                                                                           |                                                                                       |
| Siemens                                                                             | Siemens Medical:Sirecust 700, Sirecust 630                                                                                                                                                                                                                                                                                                            |             | EX0001A          | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>6p6c jack,   | 257  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0001B          | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>Din 2.0 4pin   | 257  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0001C          | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>5pin,        | 257  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0001D          | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>7p jack,         | 257  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0001E          | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>5pin          | 257  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                                                                       |             | EX0001F          | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>Round 4p    | 257  |  |
|  |                                                                                                                                                                                                                                                                                                                                                       |             |                  |                                                                        |                                                                                           |                                                                                       |



## Invasive Blood pressure Cables and transducers

| Compatible Brand          | Compatible Models No.                                                                                                             | OEM P/N | Evotec Part number | Description<br>(Name,Length,Plug)                                      | INSTRUMENT CONNECTOR                                                                       | Yoke<br>( transducers side)                                                           |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------|---------|--------------------|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Siemens<br/>Drager</b> | Siemens/Darger:SC8002XL,7000,SC8000,SC9000XL<br> |         | EX0002A            | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>6p6c jack,  | 258<br> |    |
|                           |                                                                                                                                   |         | EX0002B            | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>Din 2.0 4pin  | 258<br> |    |
|                           |                                                                                                                                   |         | EX0002C            | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>5pin,       | 258<br> |    |
|                           |                                                                                                                                   |         | EX0002D            | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>7p jack,        | 258<br> |    |
|                           |                                                                                                                                   |         | EX0002E            | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>5pin         | 258<br> |    |
|                           |                                                                                                                                   |         | EX0002F            | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 7pin>Round 4p   | 258<br> |    |
| <b>Siemens</b>            | Siemens/Darge<br>                              |         | EX00031A           | Abbott Transducer Adapter Cable, 12 ft. (3.6M), DB15M>6p6c jack,       |         |    |
|                           |                                                                                                                                   |         | EX00031B           | Utah Transducer Adapter Cable, 12 ft. (3.6M), DB15M>Din 2.0 4pin       |       |  |
|                           |                                                                                                                                   |         | EX00031C           | Edward Transducer Adapter Cable, 12 ft. (3.6M), DB15M>5pin,            |       |  |
|                           |                                                                                                                                   |         | EX00031D           | BD Transducer Adapter Cable, 12 ft. (3.6M), DB15M>7p jack,             |       |  |
|                           |                                                                                                                                   |         | EX00031E           | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), DB15M>5pin              |       |  |
|                           |                                                                                                                                   |         | EX00031F           | Philips Transducer Adapter Cable, 12 ft. (3.6M), DB15M>Round 4p        |       |  |
| <b>Drager</b>             | Drager PM8014 & PM8060<br>                     |         | EX00032A           | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>6p6c jack, |       |  |
|                           |                                                                                                                                   |         | EX00032B           | Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>Din 2.0 4pin |       |  |
|                           |                                                                                                                                   |         | EX00032C           | Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>5pin,      |       |  |
|                           |                                                                                                                                   |         | EX00032D           | BD Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>7p jack,       |       |  |
|                           |                                                                                                                                   |         | EX00032E           | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>5pin        |       |  |
|                           |                                                                                                                                   |         | EX00032F           | Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 10pin>Round 4p  |       |  |

## Invasive Blood pressure Cables and transducers

| Compatible Brand | Compatible Models No. | OEM P/N | Evotec Part number | Description<br>(Name,Length,Plug)                                         | INSTRUMENT CONNECTOR                                                                | Yoke<br>( transducers side)                                                         |
|------------------|-----------------------|---------|--------------------|---------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>Choice:</b>   |                       |         | EX0033A            | Abbott Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>6p6c jack, |  |  |
|                  |                       |         | EX0033B            | Utah Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>Din 2.0 4P   |  |  |
|                  |                       |         | EX0033C            | Edward Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>5pin,      |  |  |
|                  |                       |         | EX0033D            | BD Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>7p jack,       |  |  |
|                  |                       |         | EX0033E            | MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>6pin        |  |  |
|                  |                       |         | EX0033F            | Philips Transducer Adapter Cable, 12 ft. (3.6M), Redel 1-Rib 6p>Round 4p  |  |  |

### ADAPTER (Exchange pressure transducers type)

| Compatible Brand | Compatible Models No.                              | OEM P/N | Evotec Part number | Description<br>(Name,Length,Function)           | Plug                                                                                  | Yoke<br>( transducers side)                                                           |
|------------------|----------------------------------------------------|---------|--------------------|-------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>BD</b>        | Compatible: All monitors that accept BD plugs      |         | EX0034A            | Pressure transducers Adapter, BD to Abbott      |  |  |
|                  |                                                    |         | EX0034B            | Pressure transducers Adapter, BD to Utah        |  |  |
| <b>pvb</b>       | Compatible: All monitors that accept pvb plugs     |         | EX0035A            | Pressure transducers Adapter, pvb to Abbott     |  |  |
|                  |                                                    |         | EX0035B            | Pressure transducers Adapter, pvb to Utah       |  |  |
| <b>Philips</b>   | Compatible: All monitors that accept Philips plugs |         | EX0036A            | Pressure transducers Adapter, Philips to Abbott |  |  |
|                  |                                                    |         | EX0036B            | Pressure transducers Adapter, Philips to Utah   |  |  |



## Invasive Blood pressure Cables and transducers

### ADAPTER CABLE

| Compatible Brand    | Compatible Models No.                                                                                                         | OEM P/N     | Evotec Part number | Description (Name,Length,Plug)                                                              | INSTRUMENT CONNECTOR                                                                     | Yoke (transducers side)                                                                  |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------|-------------|--------------------|---------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| <b>Mennen</b>       |                                              |             | EX0025A            | Adapter cable, pressure transducer interface cable, 35 cm, Round 8p>AMP 6J                  | 015   | 270   |
|                     |                                              | 800-030-270 | EX0003A            | Adapter cable, pressure transducer interface cable, 35 cm, Round 10p>2* 6J                  | 260   | 269   |
| <b>Siemens</b>      |                                              | 13646       | EX0004A            | Adapter cable, pressure transducer interface cable, 35 cm, Rec16p> 2*Rou 10p                | 261   | 272   |
| <b>Spacelabs</b>    | For use with the 90496 and 690-0021-00.<br> | 700-0028-00 | EX0024A            | Adapter cable, pressure transducer interface cable, 122 cm/4 ft                             | 082   | 271   |
| <b>GE-marquette</b> | GE Solar, Dash, Eagle, Tram<br>             |             | EX0014H            | Adapter cable, pressure transducer interface cable, • Durable cable for extended cable life | 413  | 298  |

### DISPOSABLE PRESSURE TRANSDUCER

| Compatible Brand      | Product Guide                                                                       | Transducer Connector                                                                | Evotec Part number | Description                                                                                                                               |
|-----------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Abbott Medex</b>   |  |  | EX0020A            | Pressure transducer kit for patientmount application, single-use<br>Connector: RJ11-6P6C (No.: 262)<br>Requires IBP adapter cable.        |
| <b>UTAH</b>           |                                                                                     |  | EX0020B            | Pressure transducer kit for patientmount application, single-use<br>Connector: 1.58 Jack x 4pin (No.: 263)<br>Requires IBP adapter cable. |
| <b>Edwards Baxter</b> |                                                                                     |  | EX0020C            | Pressure transducer kit for patientmount application, single-use<br>Connector: 5pin (No.: 264)<br>Requires IBP adapter cable.             |
| <b>BD</b>             |                                                                                     |  | EX0020D            | Pressure transducer kit for patientmount application, single-use<br>Connector: Round 7pin (No.: 265)<br>Requires IBP adapter cable.       |
| <b>MEDEX</b>          |                                                                                     |  | EX0020E            | Pressure transducer kit for patientmount application, single-use<br>Connector: 5pin (No.: 267)<br>Requires IBP adapter cable.             |
| <b>Philips</b>        |                                                                                     |  | EX0020F            | Pressure transducer kit for patientmount application, single-use<br>Connector: 1.0 x 4pin (No.: 268)<br>Requires IBP adapter cable.       |
| <b>pvb</b>            |                                                                                     |  | EX0020G            | Pressure transducer kit for patientmount application, single-use<br>Connector: 1.5 Jack x 5pin<br>Requires IBP adapter cable.             |