


NiBp Connectors



Type Disignation

E Y □ □ □ □ □ □

Yoke code & cuff type code: T1:Thigh,Single-tube; T2:Thigh,two-tube; L1: Large Adult,Single-tube; L2: Large Adult,two-tube; A1: adult,Single-tube; A2: adult,two-tube; C1: Child,Single-tube; C2: Child,two-tube; I1: Infant,Single-tube; I2: Infant,two-tube; N1: New Bo,Single-tube; N2: New Bo,two-tube;

Instrument Connector code. If the product is cuff, the lattice is useless, but use the "000" to express. As Y000A1

Y: Indicated Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Compatibility	Product Guide	OEM P/N	Evotec Part Number	Description	Perfect fit for
Unfitted			EHL0024B	Single Tubing, Black , Inner diameter: 5/32 in. Out diameter: 6mm, Purchasing Unit: M - Material: Silicone	all Barbed diameter 5/32 in. air connector
Unfitted		GE P/N:300 959	ESL0008A	Single Tubing, Gray , Inner diameter: 5/32 in. Out diameter: 6mm, Purchasing Unit: M - Material: TPU	all Barbed diameter 5/32 in. air connector
Unfitted			ESL0008B	Two-Tubing, Gray, Inner diameter: 5/32 in. Out diameter: 6mm*2, Purchasing Unit: M - Material: TPU	all Barbed diameter 5/32 in. air connector
Air Hose Connector		GE P/N:300 670	EA 01	5/32 in. ID Barbed to 5/32 in. ID Barbed, Plastic, Packaging: 10pcs/pk	all Barbed diameter 5/32 in. air connector
L&T, Critikon, Datex (Neonatal)		Welch-allyn P/N:508 2-176	EA 26	3/32 in. ID Barbed to Male Luer Slip type, Plastic, Packaging: 10pcs/pk	Female Slip Luer (A03,A15,A16)
L&T, Critikon, Datex		GE P/N: 300667 Welch-allyn P/N: 5082-167	EA 25	5/32 in. ID Barbed to Male Slip Luer, Plastic, Packaging: 10pcs/pk	Female Slip Luer (A03,A15,A16)
L & T		GE P/N: 300667 Welch-allyn P/N: 5082-167	EA 02	5/32 in. ID Barbed to Male Slip Luer, Metal, Packaging: 10pcs/pk	Female Slip Luer (A03,A15,A16)
L & T		GE P/N: 300669 Welch-allyn P/N: 5082-169	EA 03	5/32 in. ID Barbed to Female Slip Luer, Metal, Packaging: 10pcs/pk	Male Slip Luer, (A02,A25,A26)
Philips (HP), Siemens, Datascope and Colin		GE P/N: 330059 Welch-allyn P/N: 5082-184	EA 04	Male Bayonet to 5/32 in. ID Barbed, Metal, Packaging: 10pcs/pk	Female Bayonet
		Welch-allyn P/N: 5082-184	EA 05	Male Bayonet to 5/32 in. ID Barbed, Plastic, Packaging: 10pcs/pk	Female Bayonet
		GE P/N:330 060	EA 06	5/32 in. ID Barbed to Female Bayonet, Metal, Packaging: 1piece/pk	Male Bayonet

Compatibility	Productor Guide	OEM P/N	Evotec # No.	Description	Perfect fit for
Spacelabs, Datascope and Colin		GE P/N: 3006 68 Welch-alllyn P/N: 5082 -168	EA 07	(Plastic Luer Lock) Female Locking Luer to 5/32 in. ID Barbed, Plastic, Packaging: 10pcs/pk	Male Locking Luer
		GE P/N: 3006 65 Welch-alllyn P/N: 5082 -165	EA 08	5/32 in. ID Barbed to Male Locking Luer (Plastic Luer Lock), Plastic, Packaging: 10pcs/pk	Female Locking Luer
DINAMAP		GE P/N: 3006 64 Welch-alllyn P/N: 5082 -164	EA 09	Male Screw to 5/32 in. ID Barbed, Plastic, Packaging: 10pcs/pk	Female Screw
		GE P/N: 3006 19 Welch-alllyn P/N: 5082 -161	EA 10	1/8 in. ID Barbed to Female Screw, Plastic, Packaging: 10pcs/pk	Male Screw
GE, Marquette and Protocol		GE P/N: 3300 64 Welch-alllyn P/N: 5082 -182	EA 11	(Quick Disconnect Locking type) Female Submin to 1/8 in. ID Barbed, Plastic, Packaging: 10pcs/pk	Male Submin
GE, Marquette, Datex-Ohmeda		GE P/N: 3300 90	EA 12	(Quick Disconnect Locking type) Male Submin to 5/32 in. ID Barbed, Plastic, Packaging: 10pcs/pk	Female Submin
Marquette one-tube systems		GE P/N: 3300 57 Welch-alllyn P/N: 5082 -181	EA 13	5/32 in. ID Barbed to Female Marquette locking type Connector, Plastic, Packaging: 1piece/pk	male Marquette locking type Connector
			EA 14	5/32 in. ID Barbed to male Marquette locking type Connector, Plastic, Packaging: 1piece/pk	Female Marquette locking type Connector
ADAPTER			EA 15	5/32 in. ID Barbed X2 to Female Slip Luer,, Plastic, Packaging: 10pcs/pk	Marquette Neonatal cuff connector
			EA 21	5/32 in. ID Barbed X2 to 5/32 in. ID Barbed X2, Plastic, Packaging: 10pcs/pk	5/32 in. ID Tube
Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800/850/851			EA 16	5/32 in. ID Barbed to Female Slip Luer, Plastic, Packaging: 10pcs/pk	Marquette Neonatal cuff connector

Compatibility	Productor Guide	OEM P/N	Evotec # no.	Description	Perfect fit for
Readel			A17	5/32 in. ID Barbed to Readel air socket, Plastic, Packaging: 1piece/pk	Readel air plug
			EA 18	5/32 in. ID Barbed to Readel air plug, Plastic, Packaging: 1piece/pk	Readel air socket,
Philips			EA 19	Male Bayonet socket to 5/32 in. ID Barbed, Metal, L=29mm, Packaging: 1piece/pk	Female Bayonet
			EA 20	Male Bayonet socket to 5/32 in. ID Barbed, Metal, L=38mm, Packaging: 1piece/pk	Female Bayonet
		GE P/N: 330 060	EA 06	5/32 in. ID Barbed to Female Bayonet, Metal, Packaging: 1piece/pk	Male Bayonet socket
Marquette		GE P/N: 330 057	EA 13	5/32 in. ID Barbed to Female Marquette locking type socket, Plastic, Packaging: 1piece/pk	male Marquette locking type Connector
			EA 14	5/32 in. ID Barbed to male Marquette locking type Connector, Plastic, Packaging: 1piece/pk	Female Marquette locking type socket
Eagle 4000, Tram 100/200/250/300/350, Corometrics 118			EA 22	5/32 in. ID Barbed*2 to Round Connector, Packaging: 1piece/pk	Perfect fit for two-tubes
Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120			EA 23	Rectangular Connector, Packaging: 1piece/pk	Perfect fit for two-tubes
DINAMAP, Pro, Compact, MPS Series			EA 24	Gray Shroud Connector, Packaging: 1piece/pk	Perfect fit for two-tubes