



One piece 10 Lead ECG cables with leads



Connector type



Lead Position Marked



12 Lead: I
Standard IEC

12 Lead: A
Standard AHA

S: Standard Snap



C: Clip (type RI, Grabber)




D: Din 3.0



B: Banana 4.0



Type Disignation




Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	Evotec Number	STANDARD	Defibrillation withstand	PATIENT END TERMINAL
Bosch Cambridge Custo-med DEGO Hellige Siemens	Bosch: 603 D, 603 ER, 603 S, 606 Cambridge: CM 6000, CM 6400 Custo-med: Custo-card DEGO: Cardiomed 3 Hellige: EK 36, EK 41, EK 43, EK 53, EK 56, EK 403, EK 413, EK 512, EK 512 E, EK 512 P, Cardiosmart, CardioSys, Coma Siemens: Cardiotat 1 / 3 T / 703 / 710, Cardina C, Cardirex 3T/4/6T/62, Cathcor, Megacard Elema A/B, R, R/E, Mingodoc, Mingocard 3 / 4 / 7, Mingograph 4 / 62 / 81 Mingograph 82 / 740, Mingolog 7, Mingophon 3/7/310/410/420/710/720, Miracard 1, Sirecard440/440S/460/460S/740 Siereg, Sirecard 803, Sirecust 701	291 	EVA029SNA	AHA,	No Resistance,	Snap,
			EVA029CNA	AHA,	No Resistance,	Clip,
			EVA029DNA	AHA,	No Resistance,	Din3.0,
			EVA029BNA	AHA,	No Resistance,	Banana
			EVA029SAA	AHA,	4.7kΩ Resistance,	Snap,
			EVA029CAA	AHA,	4.7kΩ Resistance,	Clip
			EVA029DAA	AHA,	4.7kΩ Resistance,	Din3.0
			EVA029BAA	AHA,	4.7kΩ Resistance,	Banana
			EVA029SNI	IEC,	No Resistance,	Snap,
			EVA029CNI	IEC,	No Resistance,	Clip,
			EVA029CNI	IEC,	No Resistance,	Din3.0,
			EVA029BNI	IEC,	No Resistance,	Banana
			EVA029SAI	IEC,	4.7kΩ Resistance,	Snap,
			EVA029CAI	IEC,	4.7kΩ Resistance,	Clip
			EVA029DAI	IEC,	4.7kΩ Resistance,	Din3.0
			EVA029BAI	IEC,	4.7kΩ Resistance,	Banana
Von Berg (hormann) Siemens	Von Berg: Bioset 3500 / 3600 / 3700 Bioset 6000 / 8000 Siemens: Cardiotat 31 / 31 S	291 	EVA029SAA	AHA,	4.7kΩ Resistance,	Snap,
			EVA029CAA	AHA,	4.7kΩ Resistance,	Clip
			EVA029DAA	AHA,	4.7kΩ Resistance,	Din3.0
			EVA029BAA	AHA,	4.7kΩ Resistance,	Banana
			EVA029SAI	IEC,	4.7kΩ Resistance,	Snap,
			EVA029CAA	IEC,	4.7kΩ Resistance,	Clip
			EVA029DAA	IEC,	4.7kΩ Resistance,	Din3.0
			EVA029BAI	IEC,	4.7kΩ Resistance,	Banana
Hellige	Hellige: EK 51	036  	EVA011SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA011CBA	AHA,	10 kΩ Resistance,	Clip
			EVA011DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA011BBA	AHA,	10 kΩ Resistance,	Banana
			EVA011SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA011CBI	IEC,	10 kΩ Resistance,	Clip
			EVA011DBI	IEC,	10 kΩ Resistance,	Din3.0
			EVA011BBI	IEC,	10 kΩ Resistance,	Banana
Avionics (Del Mar) Custo-med Del mar Ergoline Esaote Biomed Fidelity Innomed Kontron Reynolds Schiller Trentina (Ellettronica) Von Berg Bionet Welch allyn,	Avionics (Del Mar): 910 / 920 / 930 / 940 / 942 / 950 / 960 / 970 / 980 Custo-med: Custo-norm Del mar: 910,920,930,940,942,950, 960,970,980 Ergoline: Autoruler, Autoscript 6/12 Cardiette, EK 3003 / 3012, Ergoline Esaote Biomedica: EKG P80, 120 Fidelity: LP 3000 Innomed: Heart Mirror (3-Kanal) Kontron: EKG P80 Reynolds: Cardiotrack 12 Schiller: AT 1, AT 2, AT 3, AT 5, AT 6, AT10, AT60, AT104 PC, CS100, CS6, CS 6/12, AT4, AT104, AT101, AT102 Trentina: All Excel-Types Von Berg: Bioset 9000	036  	EVA021SNA	AHA,	No Resistance,	Snap,
			EVA021CNA	AHA,	No Resistance,	Clip,
			EVA021DNA	AHA,	No Resistance,	Din3.0,
			EVA021BNA	AHA,	No Resistance,	Banana
			EVA021SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA021CBA	AHA,	10 kΩ Resistance,	Clip
			EVA021DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA021BBA	AHA,	10 kΩ Resistance,	Banana
			EVA021SNI	IEC,	No Resistance,	Snap,
			EVA021CNI	IEC,	No Resistance,	Clip,
			EVA021DNI	IEC,	No Resistance,	Din3.0,
			EVA021BNI	IEC,	No Resistance,	Banana
			EVA021SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA021CBI	IEC,	10 kΩ Resistance,	Clip
			EVA021DBI	IEC,	10 kΩ Resistance,	Din3.0
			EVA021BBI	IEC,	10 kΩ Resistance,	Banana
			EVA021BBI			

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	Evotec Number	STANDARD	Defibrillation withstand	PATIENT END TERMINAL
Bosch / dimeq Cambridge Cardioline/Remco Cardionics Hellige Honeywell/Philips: Polymed Schwarzer SECA	Bosch/dimeq: 501A,501NV,503,601 503B,503D,503NV,503V,506,506D, Cambridge: VS 6, SC 1000 Cardioline: E1,1 R/3,DELTA 3/300/360 Cardionics: all Models Hellige: Corina. EK 51 Honeywell/Philips: Cardioplan, CR 300, CR 330, ER 330 Polymed: Alpha 600 / 800 Schwarzer: VS 6 SECA: 3000	038  	EVA001SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA001CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA001DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA001BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA001SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA001CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA001DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA001BAI	IEC,	4.7k Ω Resistance,	Banana
Burdick Quinton ATRIA	BURDICK (QUINTON): Eclipse 4, Eclipse 400,Eclipse 8, Eclipse 800, Eclipse 850, Eclipse 850i, Eclipse Plus, Eclipse LE II BURDICK with serial number starting with 1067 or 92304 ATRIA: ATRIA PLUS, ATRIA3000	162  	EVA005SBA	AHA,	10 k Ω Resistance,	Snap,
			EVA005CBA	AHA,	10 k Ω Resistance,	Clip
			EVA005DBA	AHA,	10 k Ω Resistance,	Din3.0
			EVA005BBA	AHA,	10 k Ω Resistance,	Banana
			EVA005SBI	IEC,	10 k Ω Resistance,	Snap,
			EVA005CBI	IEC,	10 k Ω Resistance,	Clip
			EVA005DBI	IEC,	10 k Ω Resistance,	Din3.0
			EVA005BBI	IEC,	10 k Ω Resistance,	Banana
Burdick	EK-10	039  	EVA002SNA	AHA,	No Resistance,	Snap,
			EVA002CNA	AHA,	No Resistance,	Clip,
			EVA002DNA	AHA,	No Resistance,	Din3.0,
			EVA002BNA	AHA,	No Resistance,	Banana
			EVA002SNI	IEC,	No Resistance,	Snap,
			EVA002CNI	IEC,	No Resistance,	Clip,
			EVA002DNI	IEC,	No Resistance,	Din3.0,
			EVA002BNI	IEC,	No Resistance,	Banana
Burdick	Burdick: EK-10, Elite, Elite II, E-350, E350i, E550, E560, (1st Generation Eclipse LE/LEi, with serial number starting with 1051). Burdick #007705, Stress cable with detachable leads 007706	038  	EVA001SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA001CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA001DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA001BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA001SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA001CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA001DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA001BAI	IEC,	4.7k Ω Resistance,	Banana
Fukuda Denshi	Fukuda Denshi: Alpha 500 / 600, FCP 2201	038  	EVA001SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA001CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA001DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA001BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA001SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA001CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA001DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA001BAI	IEC,	4.7k Ω Resistance,	Banana
	Fukuda Denshi Autocardiner, Cardimax FX-2111 / FX-3010 / FX-4010, FCP-2155	034  	EVA017SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA017CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA017DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA017BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA017SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA017CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA017DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA017BAI	IEC,	4.7k Ω Resistance,	Banana
FUKUDA ME	KP-500	165  	EVA024SNA	AHA,	No Resistance,	Snap,
			EVA024CNA	AHA,	No Resistance,	Clip,
			EVA024DNA	AHA,	No Resistance,	Din3.0,
			EVA024BNA	AHA,	No Resistance,	Banana
			EVA024SNI	IEC,	No Resistance,	Snap,
			EVA024CNI	IEC,	No Resistance,	Clip,
			EVA024DNI	IEC,	No Resistance,	Din3.0,
			EVA024BNI	IEC,	No Resistance,	Banana

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	MED-LINK REF NO.	STANDARD	Defibrillation withstand	PATIENT END TERMINAL
FUKUDA ME	KP-500D	165 	EVA024SEA	AHA,	20k Ω +discharge,	Snap,
			EVA024CEA	AHA,	20k Ω +discharge,	Clip,
			EVA024DEA	AHA,	20k Ω +discharge,	Din3.0,
			EVA024BEA	AHA,	20k Ω +discharge,	Banana
			EVA024SEI	IEC,	20k Ω +discharge,	Snap,
			EVA024CEI	IEC,	20k Ω +discharge,	Clip,
			EVA024DEI	IEC,	20k Ω +discharge,	Din3.0,
			EVA024BEI	IEC,	20k Ω +discharge,	Banana
	FX-031	163 	EVA014SNA	AHA,	No Resistance,	Snap,
			EVA014CNA	AHA,	No Resistance,	Clip,
			EVA014DNA	AHA,	No Resistance,	Din3.0,
			EVA014BNA	AHA,	No Resistance,	Banana
			EVA014SNI	IEC,	No Resistance,	Snap,
			EVA014CNI	IEC,	No Resistance,	Clip,
			EVA014DNI	IEC,	No Resistance,	Din3.0,
			EVA014BNI	IEC,	No Resistance,	Banana
	Fukuda M-E: Cardisunny alpha 1000, Cardisunny alpha 6600 A / 6600 AX	038 	EVA001SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA001CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA001DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA001BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA001SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA001CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA001DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA001BAI	IEC,	4.7k Ω Resistance,	Banana

KANZ PC-104, Cardioline remco	Kanz: 103, 106 Cardioline/remco: delta 1,3,30,60 delta 1 plus, delta 3 plus, delta 60 plus	024 	EVA016SNA	AHA,	No Resistance,	Snap,
			EVA016CNA	AHA,	No Resistance,	Clip,
			EVA016DNA	AHA,	No Resistance,	Din3.0,
			EVA016BNA	AHA,	No Resistance,	Banana
			EVA016SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA016CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA016DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA016BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA016SNI	IEC,	No Resistance,	Snap,
			EVA016CNI	IEC,	No Resistance,	Clip,
			EVA016DNI	IEC,	No Resistance,	Din3.0,
			EVA016BNI	IEC,	No Resistance,	Banana
			EVA016SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA016CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA016DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA016BAI	IEC,	4.7k Ω Resistance,	Banana
Kanz PC-109	ECG 108/110/1203,1205	033 	EVA015SNA	AHA,	No Resistance,	Snap,
			EVA015CNA	AHA,	No Resistance,	Clip,
			EVA015DNA	AHA,	No Resistance,	Din3.0,
			EVA015BNA	AHA,	No Resistance,	Banana
			EVA015SNI	IEC,	No Resistance,	Snap,
			EVA015CNI	IEC,	No Resistance,	Clip,
			EVA015DNA	IEC,	No Resistance,	Din3.0,
			EVA015BNA	IEC,	No Resistance,	Banana
MARQUETTE (GEMS)	Hellige: MicroSmart, MAC 500 / 1100 MAC 1200 / 1200 ST	029 	EVA008SBA	AHA,	10 k Ω Resistance,	Snap,
			EVA008CBA	AHA,	10 k Ω Resistance,	Clip
			EVA008DBA	AHA,	10 k Ω Resistance,	Din3.0
			EVA008BBA	AHA,	10 k Ω Resistance,	Banana
			EVA008SBI	IEC,	10 k Ω Resistance,	Snap,
			EVA008CBI	IEC,	10 k Ω Resistance,	Clip
			EVA008DBI	IEC,	10 k Ω Resistance,	Din3.0
			EVA008BBI	IEC,	10 k Ω Resistance,	Banana

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	Evotec Number	STANDARD	Defibrillation withstand	PATIENT END TERMINAL
MEDICAL RESEARCH LABS (MRL)	MEDICAL RESEARCH LABS (MRL): Portable Intensive Care (PIC) System	223 	EVA012SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA012CBA	AHA,	10 kΩ Resistance,	Clip
			EVA012DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA012BBA	AHA,	10 kΩ Resistance,	Banana
			EVA012SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA012CBI	IEC,	10 kΩ Resistance,	Clip
			EVA012DBI	IEC,	10 kΩ Resistance,	Din3.0
MEDTRONIC PHYSIO CONTROL	MEDTRONIC PHYSIO CONTROL: Lifepak 12	224 	EVA013SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA013CBA	AHA,	10 kΩ Resistance,	Clip
			EVA013DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA013BBA	AHA,	10 kΩ Resistance,	Banana
			EVA013SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA013CBI	IEC,	10 kΩ Resistance,	Clip
			EVA013DBI	IEC,	10 kΩ Resistance,	Din3.0
Mindray	12-Ld module		EVA020SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA020CBA	AHA,	10 kΩ Resistance,	Clip
			EVA020DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA020BBA	AHA,	10 kΩ Resistance,	Banana
			EVA020SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA020CBI	IEC,	10 kΩ Resistance,	Clip
			EVA020DBI	IEC,	10 kΩ Resistance,	Din3.0
MORTARA	MORTARA INSTRUMENT:ELI 100 MORTARA INSTRUMENT:ELI 200 MORTARA INSTRUMENT:ELI 50	158 	EVA025SAA	AHA,	4.7kΩ Resistance,	Snap,
			EVA025CAA	AHA,	4.7kΩ Resistance,	Clip
			EVA025DAA	AHA,	4.7kΩ Resistance,	Din3.0
			EVA025BAA	AHA,	4.7kΩ Resistance,	Banana
			EVA025SAI	IEC,	4.7kΩ Resistance,	Snap,
			EVA025CAI	IEC,	4.7kΩ Resistance,	Clip
			EVA025DAI	IEC,	4.7kΩ Resistance,	Din3.0
MORTARA	Martara Surveyor 		EVA035SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA035CBI	IEC,	10 kΩ Resistance,	Clip,
NEC	YY-1063	027 	EVA003SBA	AHA,	10 kΩ Resistance,	Snap,
			EVA003CBA	AHA,	10 kΩ Resistance,	Clip
			EVA003DBA	AHA,	10 kΩ Resistance,	Din3.0
			EVA003BBA	AHA,	10 kΩ Resistance,	Banana
			EVA003SBI	IEC,	10 kΩ Resistance,	Snap,
			EVA003CBI	IEC,	10 kΩ Resistance,	Clip
			EVA003DBI	IEC,	10 kΩ Resistance,	Din3.0
Nihon Konden Biocare邦健 DongJiang东江	Nihon Konden: ECG-9620/9010/9020 ECG-9022/9110/9130/9132 ECG-8420(Except 9320) BIOCARE邦健: ECG-101/101G/300 ECG-300G,ECG-1200, ECG-9803. CMICS: ECG-11A, ECG-11B,ECG-11D, ECG-32A, ECG-32A(II), ECG-120 ,	038 	EVA001SNA	AHA,	No Resistance,	Snap,
			EVA001CNA	AHA,	No Resistance,	Clip,
			EVA001DNA	AHA,	No Resistance,	Din3.0,
			EVA001BNA	AHA,	No Resistance,	Banana
			EVA001SNI	IEC,	No Resistance,	Snap,
			EVA001CNI	IEC,	No Resistance,	Clip,
			EVA001DNI	IEC,	No Resistance,	Din3.0,
Nihon Kohden	BSM-2301, BSM-2353, BSM-5100,		EVA001BNI	IEC,	No Resistance,	Banana
			EVA036SNI	IEC,	No Resistance,	Snap,
			EVA036CNI	IEC,	No Resistance,	Clip,
			EVA036DNI	IEC,	No Resistance,	Din3.0,
			EVA036BNI	IEC,	No Resistance,	Banana

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	Evotec Number	STANDARD	Defibrillation withstand	PATIENT END TERMINAL
Nihon Kohden	Nihon Kohden: Cardiofax6151/6353, Cardiofax EKG-8270/8350/8370/8420 8423/8450/8453/8820/8823/8830/9130 Cardiofax-Q	038 	EVA001SAA	AHA,	4.7k Ω Resistance,	Snap,
			EVA001CAA	AHA,	4.7k Ω Resistance,	Clip
			EVA001DAA	AHA,	4.7k Ω Resistance,	Din3.0
			EVA001BAA	AHA,	4.7k Ω Resistance,	Banana
			EVA001SAI	IEC,	4.7k Ω Resistance,	Snap,
			EVA001CAI	IEC,	4.7k Ω Resistance,	Clip
			EVA001DAI	IEC,	4.7k Ω Resistance,	Din3.0
			EVA001BAI	IEC,	4.7k Ω Resistance,	Banana
	Cardiofax Q ECG-9110K Cardiofax Q ECG-9130K ECG-9130P, ECG-9620P special screw, imperial standard(SAE)	038 	EVA009SNA	AHA,	No Resistance,	Snap,
			EVA009CNA	AHA,	No Resistance,	Clip,
			EVA009DNA	AHA,	No Resistance,	Din3.0,
			EVA009BNA	AHA,	No Resistance,	Banana
			EVA009SNI	IEC,	No Resistance,	Snap,
			EVA009CNI	IEC,	No Resistance,	Clip,
			EVA009DNI	IEC,	No Resistance,	Din3.0,
			EVA009BNI	IEC,	No Resistance,	Banana
	ECG-9320 ECG-9522P (Original part No.: BR-911D)	374 	EVA010SCA	AHA,	20 k Ω Resistance,	Snap,
			EVA010CCA	AHA,	20 k Ω Resistance,	Clip
			EVA010DCA	AHA,	20 k Ω Resistance,	Din3.0
			EVA010BCA	AHA,	20 k Ω Resistance,	Banana
			EVA010SCI	IEC,	20 k Ω Resistance,	Snap,
			EVA010CCI	IEC,	20 k Ω Resistance,	Clip
			EVA010DCI	IEC,	20 k Ω Resistance,	Din3.0
			EVA010BCI	IEC,	20 k Ω Resistance,	Banana
PETAS		222 	EVA023SNA	AHA,	No Resistance,	Snap,
			EVA023CNA	AHA,	No Resistance,	Clip,
			EVA023DNA	AHA,	No Resistance,	Din3.0,
			EVA023BNA	AHA,	No Resistance,	Banana
			EVA023SNI	IEC,	No Resistance,	Snap,
			EVA023CNI	IEC,	No Resistance,	Clip,
			EVA023DNI	IEC,	No Resistance,	Din3.0,
			EVA023BNI	IEC,	No Resistance,	Banana
PHILIPS (HP)	PHILIPS MEDICAL SYSTEMS: M1770A PageWriter 200i , M1771A PageWriter 200 M1772A PageWriter 100	035 	EVA018SCA	AHA,	20 k Ω Resistance,	Snap,
			EVA018CCA	AHA,	20 k Ω Resistance,	Clip
			EVA018DCA	AHA,	20 k Ω Resistance,	Din3.0
			EVA018BCA	AHA,	20 k Ω Resistance,	Banana
			EVA018SCI	IEC,	20 k Ω Resistance,	Snap,
			EVA018CCI	IEC,	20 k Ω Resistance,	Clip
			EVA018DCI	IEC,	20 k Ω Resistance,	Din3.0
			EVA018BCI	IEC,	20 k Ω Resistance,	Banana
Spacelabs	90496 Ultraview Module	022 	EVA027SBA	AHA,	10 k Ω Resistance,	Snap,
			EVA027CBA	AHA,	10 k Ω Resistance,	Clip
			EVA027DBA	AHA,	10 k Ω Resistance,	Din3.0
			EVA027BBA	AHA,	10 k Ω Resistance,	Banana
			EVA027SBI	IEC,	10 k Ω Resistance,	Snap,
			EVA027CBI	IEC,	10 k Ω Resistance,	Clip
			EVA027DBI	IEC,	10 k Ω Resistance,	Din3.0
			EVA027BBI	IEC,	10 k Ω Resistance,	Banana
Shanghai Konden 上海光电	ShangHai Kohden: ECG-6511, 6151, 6551, 8110P/K, 6353, 6403, 6504	028 	EVA004SNA	AHA,	No Resistance,	Snap,
			EVA004CNA	AHA,	No Resistance,	Clip,
			EVA004DNA	AHA,	No Resistance,	Din3.0,
			EVA004BNA	AHA,	No Resistance,	Banana
			EVA004SNI	IEC,	No Resistance,	Snap,
			EVA004CNI	IEC,	No Resistance,	Clip,
			EVA004DNI	IEC,	No Resistance,	Din3.0,
			EVA004BNI	IEC,	No Resistance,	Banana
ZhonLian JianYi 中联健益	JIANYI: XD-106,XD-104B/D		EVA031SNI	IEC,	No Resistance,	Snap,
			EVA031CNI	IEC,	No Resistance,	Clip,
			EVA031DNI	IEC,	No Resistance,	Din3.0,
Trentina:	All Excel-Types 	294 	EVA031BNI	IEC,	No Resistance,	Banana
ZOLL			EVA031BBI	IEC,	10 k Ω Resistance,	Banana