

Electrosurgical accessories


P □□□-□□□-□□

Length code, as 13 indicated 13 ft. ; 09 indicated 9 ft

Patient end type code.

Connector code: 143 indicated Din 1.5 jack;, 134 indicated Din 2.0 plug, 112 indicated 3.5 mono plug

P: Indicated Product Categories: ESU pencil and Patient Return Plate Cable

MED-LINK REF NO.	Productor Guide	Description, Name, Length, Connector type
EP281-280-13	 280 281	Patient Return Plate Cable, Length: 13 ft. (4 m) PLUG: 2.35 jack
EP053-281-13	 280 282	Patient Return Plate Cable, Length: 13 ft. (4 m) PLUG: 6.35 Mono Plug
EP053-296-10	 296 282	Patient Return Plate Cable, Length: 10 ft. (3 m) PLUG: 6.35 Mono Plug to Banana plug
EP283-284-10	 284 283	Bipolar Electrode Cable, Length: 10 ft. (3 m) PLUG: 2*Banana Plug Fit Euro style bipolar electrodes
EP351-148-13	 148 351	Connecting cable for electrosurgery Original No. : 20192-117 Intened use : the connecting cables are intended for the connection of accessories to electrosurgical equipment. Connector type: 7.9 mono plug to 3.0 Jack. Length: 13feet; Clour: Grey
EP296-297-10	 297 296	Bipolar Electrode Cable, Length: 10 ft. (3 m) PLUG: 2*Banana Plug, Fit Euro style bipolar electrodes
EP296-315-10	 315 296	Bipolar Electrode Cable, Length: 10 ft. (3 m) PLUG: 2*Banana Plug, Jack: $\phi 2.0 \times 2$ Fit USA style bipolar electrodes
EP285-286-09	 285 286	product name: Disposable Electrosurgical pencil, Cautery pencil;Categories: Surgical apparatus;Length of handle(mm):170mm; Length of blade(mm):40mm; Diameter of blade:2.3mm; Length of wire(m): 9 ft (2.6m); Package way:100pcs/box; Size of box: 45*35*35cm
EP285-286-09-R	 285 286	product name: Reusable Electrosurgical pencil, Cautery pencil;Categories: Surgical apparatus;Length of handle(mm):170mm; Length of blade(mm):40mm; Diameter of blade:2.3mm; Jacket of Cable: silicone, Length of wire(m): 9 ft (2.6m); Package way:1piece/pk;
EP287		Disposable Electrosurgical Patient Plate-Monopolar, Size: 200x105mm
EP288		Disposable Electrosurgical Patient Plate-Bipolar , Size: 200x105mm