







Incubator accessories

Line #	Part Number	Description	Rand Excl VAT	Picture
	FILTERS			
	EVO-Atom - S	ATOM transport incubator Filter 195mm X 90 mm		
	EVO - Atom-L	ATOM V-80, V-85, V850 series 165mm X 510mm		
	ET-94570 Replaces Air Shields #26-945-70	For use with: Airshield Incubators C-86,100,200,450,500,550 620mm x 150mm		
	ET-10130 Replaces Air Shields #83-101-30	Airshield Isolette C-2000 series 150mm x 80 mm		
	ET-94570	Ningbo David Incubators 620mm x 150mm		
	EVO-TI2000	Ningbo David Transport incubator T2000 190mm x 70mm		
	EVO-Medix	MEDIX incubators 850mm X 1000mm		
	SLEEVES & CUFFS			

	ET-93481	Air shields ,ATOM, Ningbo David , Dison, High strength - IRIS		
	ET-03300	Ohmeda/Ohio Models 190A, 190 ASC and Care Plus Systems High strength - IRIS		
	ET-61500	Air shields ,ATOM, Ningbo David , Dison, High strength - OVAL		
	ET-80217	Color: Assorted Pink & Blue		
	MATTRESS			
	Standard	Mattress - Ningbo , Airshields , ATOM , Standard size 610 X 330 X 30mm CUSTOM SIZES AVAILABLE		
	Medix	Mattress for Medix TI		
	REUSABLE Temperature sensors			
	ET-10241 Replaces Air Shields #83-102-41 (Drager # MU 12533 Preemie Probe 5)	Air Shields C2000 Isolette Incubators and Versalet 7700 Care Centers.		
	ET-08750 Replaces Ohmeda 6600-0875-700	All Ohmeda Giraffe Incubators.		
	ET-11225 Replaces Air Shields #81-300-15 (Drager # MU 11225 Preemie Probe 3)	Air Shields Resuscitaire Infant Radiant Warmers.		

	ET-20970 Replaces Air Shields #68-209-70	Air Shields Incubators C-100, C-100QT, C-450 QT-550 Infant Intensive Care Systems: 78-000-50, 78-001-70, 78-280-70, 78-070-70 Infant Warmers: 78-270-70, 78-275-70 Delivery Room Warmers: 78-285-75, 78-285-70 Birthing Room Warmer: 78-290-70, and IICS-90 Series.		
	ET-20980 Replaces Air Shields #68-209-80	Air Shields C-100 & C-200 Incubators.		
	ET-57514 Replaces Air Shields #17-575-14	Air Shields models: C-86 Incubator and SC-78 Infant Warmer.		
	ET-009 Replaces Air Shields #BCZ009	Air Shields Stabilet models 1250, 1500, 2000C and 2000D Infant Radiant Warmers.		
	ET-97700 Replaces Ohmeda/Ohio #208-0697-700, LA003 and LA005	Incubators: Ohmeda/Ohio IC, GC Care Plus. Warmers: Ohmeda Panda Baby Warmers and Ohmeda 3000, 4000 & 5000 Series of Infant Warmer Systems.		
	ET-54700 Replaces Ohmeda/Ohio #208-0754-700 and LA 001	Ohmeda/Ohio 190, 190A and 190ASC Incubators & Infant Warmers.		
	Disposable Temp Probes			
	ET-08730 Replaces Ohmeda 6600-0873-700 and 6600-0874-700	All Ohmeda Giraffe Incubators and the new Panda Warmer.		
	ET-10415 Replaces Air Shields #83-104-15 (Drager MU 12551 Preemie Probe 5-L)	Air Shields C2000 Isolette Incubators and Versalet 7700 Care Centers. This is a long-wire disposable probe, no other wires needed.		

	ET-11000 Replaces Draeger MX 11000	Draeger Caleo systems		
	ET-100 Replaces Air Shields Preemie Probe I, #58-209-20 and #68-209-30			
	ET-86 Replaces Air Shields Preemie Probe II, #68-209-25 and #68-209-35			
	ET-3000 Replaces Ohmeda/Ohio #208-0697-700, LA003, and LA005			
	Various parts			
	E- Grom	IV grommet – fit most incubators		
	E-Schock	Shock absorber for incubator motors		
	E-Bearing	Motor bearings for Airshields motors		
	E-Fan	Incubator motor fan (while stocks last)		

	Inc. Motor-A	Incubator motor for Airshields		
	Inc. Motor-N	Incubator motor for Ningbo David		
	Photo-20-52	Phototherapy Lights - BLUE		
	Photo-20-54	Phototherapy Lights – White		
	BIL- Test	Jaundice light meter		
		WE ALSO SUPPLY VARIOUS OTHER INCUBATOR PARTS – NOT LISTED HERE		