


NiBp Tubing with connectors



Type Disignation



















E Y □ □ □ □ □ □

Yoke code & cuff type code: T1:Thigh,Single-tube; T2:Thigh,two-tube; L1: Large Adult,Single-tube; L2: Large Adult,two-tube; A1: adult,Single-tube; A2: adult,two-tube; C1: Child,Single-tube; C2: Child,two-tube; I1: Infant,Single-tube; I2: Infant,two-tube; N1: New Bo,Single-tube; N2: New Bo,two-tube;

Instrument Connector code. If the product is cuff, the lattice is useless, but use the "000" to express. As Y000A1

Y: Indicated Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Compatible Brand	Compatible Models No.	OEM P/N	Evotec Part No.	Description (Name,Length,Plug)	Instrument Connector	Yoke (cuff side)
Philips	Connects all sizes of Philips' reusable and disposable adult and pediatric NIBP cuffs to the monitor. Use with adult and pediatric cuffs only. CANNOT BE USED WITH NEONATAL CUFFS. Use with all Philips models except A1/A3	M1598B	EYA25A06-05	Pressure Cuff Interconnect Tubing, Single tube, 5ft. 1.5M, 6.0 Air plug > Bayonet type, perfect fit for 5082-184 cuff connector	347 	348 
		M1599B	EYA25A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10ft. 3.0M, 6.0 Air plug > Bayonet type, perfect fit for 5082-184 cuff connector	347 	348 
Philips	NIBP hose connects to 40401X, M45xx, M157x series pressure cuffs. Use with Philips Instrument: A1/A3, M3921A, M3923A, M3926A, M3927A, M3928A, M3929A	M3918A	EYA06A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10 ft 3.0 M, Bayonet type > Bayonet type perfect fit for 5082-184 cuff connector	348 	348 
Simens Mindray Spacelabs MEK,NPB , Contec	Spacelabs:90496 Module, Contec: CMS9000	Spacelabs P/N: 714-0018-00	EYA06A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10 ft 3.0 M, Bayonet type > Bayonet type perfect fit for 5082-184 cuff connector	348 	348 
Goldway			EYA10A10-12	Blood Pressure Tubing (Air Hoses), single-tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,		
Creative科瑞康 Neusoft东软 Jeteem井田医电	Jeteem:BD6000		EYA18A01-10	Pressure Cuff Interconnect Tubing, Single tube, 3.0M, Readel Air plug> Barbed to barbed connector	319 	
Colin, Criticare, CSI Datascope, NihonKohden, Datascope			EYA07A08-10	Pressure Cuff Interconnect Tubing, Single tube, 3.0M, Male Plastic Luer Lock> Female Plastic Luer Lock		
GE Dinamap, Pro, MPS	DINAMAP, Pro, Compact, MPS Series	107365	EYA24A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray Shroud to Screw Connector, Gray		
		107366	EYA24A10-24	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 24 ft., Gray Shroud to Screw Connector,		
		107368	EYA24A16-12	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 12 ft., Gray Shroud to Female Slip Luer,		
		107363	EYA24A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray Shroud to Submin Connector, Gray		

Compatible Brand	Compatible Models No.	OEM P/N	Evotec Part No.	Description (Name,Length,Plug)	Instrument Connector	Yoke (cuff side)
GE marquette	Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800, 850/851, Corometrics 120	9461-217	EYA23A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft. Gray, Rectangular Connector to Screw Connector,		
		414873-001	EYA23A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray, Rectangular Connector to Submin Connector,		
		414874-001	EYA23A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Rectangular Connector to Female Slip Luer,		
GE medical Critikon, Kontron, Invivo	DINAMAP 1846, 8100, 9300 and Plus Series Monitors	8841	EYA10A10/2-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,		
		88845	EYA10A10/2-24	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 24 ft., Screw Connector to Screw Connector,		
		8841	EYA10A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Screw Connector to Female Slip Luer,		
		88845	EYA10A12-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Screw Connector to Female Slip Luer,		
GE medical Corometrics GE Datex	Eagle 4000, Tram 100/200/250/300/350, Corometrics 120	9461-203	EYA22A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Submin Connector,		
		9461-214	EYA22A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,		
		9461-205	EYA22A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Round Connector to Female Slip Luer	