


Patient Cables Two piece (Trunk + Fly leads)



Instrument Connector



Yoke type and code

Din 3-Ld, Code: C	Din 5-Ld, Code: C	Colin, 5-Ld, Code: M
Lifelinc(LL), 3-Ld, Code: B	Lifelinc(LL), 5-Ld, Code: B	TRU-link, 5-Ld, Code: L
ASP, 3-Ld, Code: H	ASP, 5-Ld, Code: H	VS, 5-Ld, Code: E
AAMI(AA), 3-Ld, Code: D	AAMI(AA), 5-Ld, Code: D	NK, 5-Ld, Code: N
EURO, 3-Ld, Code: J	EURO, 5-Ld, Code: J	NK, 6-Ld, Code: P
Conmed, 3-Ld, Code: G	Conmed, 5-Ld, Code: G	Multi-Link, 5-Ld, Code: Q
NEC, 3-Ld, Code: K	NEC, 3-Ld, Code: K	Lohmeier multi-link: JM
Multi-link, Code: R	Fukuda 16p, Code: Y	Philips M1008A yoke: DG


Type Disignation

E E □ □ □ □ □ □





Color coding: **7A**: 7-Ld, AHA(AAMI); **10A**: 10-Ld, AHA; **7I**: 7-Ld, IEC; **10I**: 10-Ld, IEC; **5A**: 5-Ld, AHA(AAMI); **3H**: 4-Ld, AHA; **3A**: 3-Ld, AHA; **5I**: 5-Ld, IEC; **3IH**: 4-Ld, IEC; **3I**: 3-Ld, IEC; **5AI**: 5-Ld, AHA/IEC; **3HIH**: 4-Ld, AHA/IEC.
 Patient end termination type: **S**: snap; **C**: Clip (grabber).
 If the product is a cable, the lattice is useless, but use the "I" to express, as ED008-3A
 Brand code, as 041 indicated Bionet; if the product is a leads, the lattice is indicated length, as 040 indicated 40 in.
 Yoke type: **A**: one-piece; **B**: Lifelinc, (LL); **C**: Din style; **D**: AA style; **E**: VS style, other sytle or details as above
 E: Indicated patient minitor ECG cable and lead wires

Evotec product

Compatible Brand	Compatible Models No.	OEM P/N	Evotec Part Number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
Din Style Normal	<u>AIR SHIELDS (NARCO, Vickers)</u> :HRRM-71-1/2 ,Porta-Fib I, II, III, System 5 (new), System 6 (Neonatal); <u>ATL ULTRASOUND</u> : 800 ,800 Plus, Apogee , <u>CX ATL ULTRASOUND</u> : CX200, CX400, HDI 9/2000, Interspec Cardiology XL, Mark IV, Mark VII , Ultramark 4 , Ultramark 4 Plus, Ultramark 5 , Ultramark 7, Ultramark 9; <u>BIOCHEM INTERNATIONAL</u> (BCI):3101, 3404, 6100, 6200, 9100, 9200, Advisor , Autocorr <u>COLIN Medical instruments</u> : BP308, BP408, CBM3000 <u>CRITICARE SYSTEMS, INC.</u> (CSI):504US, 602-4USP, Poet 2, Scholar <u>GE medical systems</u> : Corometrics 556, Critikon 8100T, Critikon Dinamap 8700 Series, Critikon Dinamap 9700 Series <u>IVY biomedical systems</u> :101 (before 1989), 402, 700	D-1340	EEC020-3A/I	Din 3Ld Patient Cable, 9ft AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires	004 	101 
		D-1540	EEC020-5A/I	Din 5Ld Patient Cable, 9ft AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires	004 	102 
			EEC021-3A/I	Din 3Ld Patient Cable, 9ft AHA/IEC, straight 6pin>Din 3J, 1KΩ resistance fit for Din style LeadWires	007 	101 
			EEC021-5A/I	Din 5Ld Patient Cable, 9ft AHA/IEC, straight 6pin>Din 5J, 1KΩ resistance fit for Din style LeadWires	007 	102 
Lifeline (LL) Style Normal	<u>LIFE SCIENCE INSTRUMENTS</u> (LSI):LS5 , 261 <u>MEDICAL DATA electronics</u> (MDE):Escort 100, Escort 101, Escort 102, Escort 200, Escort 300, Prism , 360SL, 450SL, 521A, PortaPak 80 , PortaPak 90 , PortaPak 1000; <u>Medtronic Physio-control</u> : Lifepak 6 , VSM-1ESF (w/resistors); <u>NELLCOR</u> :N-200 C-LOCK, N-250 C-LOCK, N-1000 , NPB-3900, NPB-4000; <u>NIHON KOHDEN</u> :BSM-8301A, BSM-8302A, BSM-8500A , BSM-8502A, OEC-8102A (s/n>20156), OEC-8105A (s/n>20226) , OEC-7102A , OEC-8108A, TEC-6200A (s/n>20421) , TEC-7100A , TEC-7200A , TEC-7300A; <u>PACE-TECH</u> :Minipack 300 , Minipack 3000, Minipack 3100, MiniMax 4000CL, Vitalmax 800 Plus, Vitalmax 2100, Vitalmax 2200, Vitalmax 4000, Vitalmax 4000CL, Vitalmax 4000 Modular, Vitalmax 4100, Vitalmax 4100CL;		EEB020-3A/I	LL 3Ld Patient Cable, 9ft AHA/IEC, straight 6pin>LL 3LD, no resistance fit for LL style LeadWires	004 	107 
			EEB020-5A/I	LL 5Ld Patient Cable, 9ft AHA/IEC, straight 6pin>LL 5LD, no resistance fit for LL style LeadWires	004 	105 
			EEB021-3A/I	LL 3Ld Patient Cable, 9ft AHA/IEC, straight 6pin>LL 3LD, 1KΩ resistance fit for LL style LeadWires	007 	107 
			EEB021-5A/I	LL 5Ld Patient Cable, 9ft AHA/IEC, straight 6pin>LL 5LD, 1KΩ resistance , fit for LL style LeadWires	007 	105 
Conmed Style Normal	<u>PHILIPS MEDICAL SYSTEMS</u> : M3921A A1 Compact Monitor , M3923A A1 Compact Monitor, M3925A A1 Compact Monitor , M3927A A3 Compact Monitor , M3929A A3 Compact Monitor <u>PIONEER</u> (American Optical): 26125 , 33167, 5260 , 5330 , Pulsar 4; <u>PPG</u> (Datamedix, Becton-Dickinson):Guardian 1001, Guardian N1001 <u>SPACELABS BURDICK</u> :Medio 3, Medio 4 , Medio 5 , 4045 Defib, 500 Series, 600 Series, 900 Series (w/o respiration) , 90300 Series, 90400 Series, 90600 Series , 90700 Series, TEK-208, TEK-408, TEK-511 , TEK-512 , TEK-514, TEK-521, PC Bedside , PCMS System; <u>TEKTRONIX (Vitatek)</u> :4045 Defib , 208 , 400 Series, 500 Series; <u>WELCH ALLYN</u> (Protocol):Propaq 100 Series, Propaq 200 Series, Propaq Encore, Propaq CS , Ultra Smart, Atlas <u>ZOLL</u> :D900, PD1200, PD1400		EEG020-3A/I	Conmed 3Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 3LD, no resistance, fit for Conmed style LeadWires	004 	244 
			EEG020-5A/I	Conmed 5Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 5LD, no resistance, fit for Conmed style LeadWires	004 	245 
			EEG021-3A/I	Conmed 3Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 3LD, 1KΩ resistance, fit for Conmed style LeadWires	007 	244 
			EEG021-5A/I	Conmed 5Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 5LD, 1KΩ resistance, fit for Conmed style LeadWires	007 	245 

Compatible Brand	Compatible Models No.	OEM P/N	Evotec part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
Bionet Korea	BM3		EEB041-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 8pin>LL-3LD, 1KQ resistance fit for LL style LeadWires	500 	107 
			EEB041-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 8pin>LL-5LD, 1KQ resistance fit for LL style LeadWires	015 	105 
Biosys Korea	BPM-103		EEC040-3A/I	3Ld Patient Cable,9ft AHA/IEC, Readlel 6pin>Din 3J, 1KQ resistance fit for Din style LeadWires	014 	101 
			EEB040-5A/I	5Ld Patient Cable,9ft AHA/IEC, Readlel 6pin>Din 5J, 1KQ resistance fit for Din style LeadWires	014 	102 
Colin	BP88 / BP306		EEM062-5I	5Ld Patient Cable,9ft IEC, L type round 6pin>6J, no resistance fit for Colin LeadWires	013 	083 
	BP88S		EEM061-5I	5Ld Patient Cable,9ft IEC, straight 6P>6J, no resistance fit for Colin LeadWires	025 	083 
Creative科瑞康 Spurt世博达 Biolight宝莱特	Creative: PC9000S,PC9000,UP6000, Biolight: M12		EEC046-3A/I	3Ld Patient Cable,9ft AHA/IEC, Readlel 6pin>Din 3J, 1KQ resistance fit for Din style LeadWires	014 	101 
			EEC046-5A/I	5Ld Patient Cable,9ft AHA/IEC, Readlel 6pin>Din 5J, 1KQ resistance fit for Din style LeadWires	014 	102 
Datascope	2000,2001,2002,2100,2200,3000, 3000A ,Ballon Pump,M/D3 w/lead select,M/D3A,M/D4,Passport,Passport IL,Passport XG,Passport XL,Passport 2,Passport 2LT,Spectrum System 90, System 90T , System 96, System 97		EEC002-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires	005 	101 
			EEC002-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires	005 	102 
Datex	AS/3, CS/3, S/5, Cardiocap 1, Cardiocap 2 , Cardiocap 5 , Light Monitor , Satliteplus		EEB010-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires	023 	107 
			EEB010-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires	023 	105 

Compatible Brand	Compatible Models No.	OEM P/N	Evotec part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
Datex-Ohmeda	Datex-Ohmeda Patient monitor		EEH010-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 3LD, 4.7KΩ resistance fit for ASP style LeadWires	023 	183 
			EEH010-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 5LD, no resistance fit for ASP style LeadWires	023 	221 
Drager Siemens	Drager UM1/3/3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog8000/8000 plus, PB8800,Serie 8000, Cicero, Parameterbox,		EEC057-3A/I	SG 3Ld Patient Cable,9ft AHA/IEC, Round 10pin>Din 3J, no resistance fit for Din style LeadWires	170 	101 
	Multimed-Pot system fur serie SC 6000, SC 6002XL, SC 7000, SC 9000, Art. Nr.3368391 (8.2ft) Art. Nr. 5950196 (4.9ft)	3368391	EEQ056-5A/I	SC9000XL Multi-link Cable, 8.2ft AHA/IEC, 1KΩ resistance fit for Euro style LeadWires	178 	179 
Siemens	Sirecust Series 400, 500, 600, Serie 700, 900, 1200, Sirecust 341/610/620/630/720/722, Sirecust 730/732/900/960/961/1260, Sirecust 1281/1801/400/-500, ECG Monitor, Polymed Modul, Modulbox Sirem		EEB052-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires	017 	107 
			EEB052-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires	017 	105 
			EEJ052-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>Euro 3LD, no resistance fit for Euro style LeadWires	017 	242 
			EEJ052-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>Euro 5LD, no resistance fit for Euro style LeadWires	017 	243 
GE-Medical Marquette	GEMS:Dash PRO 4000, Dash PRO 3000, Dash PRO 2000, Dash 1000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000, Solar 3000, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems		EEE007-5A	5Ld Patient Cable,9ft AHA, Round 10J>VS 5LD, no resistance fit for VS style LeadWires	020 	089 
			EEE007-5I	5Ld Patient Cable,9ft IEC, Round 10J>VS 5LD, no resistance fit for VS style LeadWires	020 	089 
GE-Hellige	Hellige,Servomed 4, Generation, SCP 851, SCP 852, SCP 844, SCS 804,, CardioServ, VICOM-sm		EEB011-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10P>LL 3LD, no resistance fit for LL style LeadWires	021 	107 
			EEB011-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10P>LL 5LD, no resistance fit for LL style LeadWires	021 	105 

OEM	OEM MODELS NUMBER	OEM P/N	Evotec part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
GE-Hellige	Hellige,Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844, SCS 804,, CardioServ , VICON-sm		EEH011-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires		 183
			EEH011-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires		 221
GE-medical Marquette	Pro1000		EEH021-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6pin>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires		 183
			EEH021-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6pin>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires		 221
Goldway	Goldway: UT4000A,UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800		EEB002-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6p>LL 3LD, no resistance fit for LL style LeadWires		 107
			EEB002-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6p>LL 5LD, no resistance fit for LL style LeadWires		 105
Kontron	7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, KAAAT, Kolormon, Micromon, Minimon, Supermon , Vitalmon		EEB031-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 12P>LL 3LD, no resistance fit for LL style LeadWires		 107
			EEB031-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 12P>LL 5LD, no resistance fit for LL style LeadWires		 105
MEK Korea	MP1000,MP600,MP500		EEC029-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, DB9M>Din 3J, no resistance fit for Din style LeadWires		 101
			EEC029-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, DB9M>Din 5J, no resistance fit for Din style LeadWires		 102
	M160246, M360697		EEC033-3A/I	3Ld Patient Cable,9ft AHA/IEC, MD6M>LL 3LD, no resistance fit for LL style LeadWires		 107
			EEB033-5A/I	5Ld Patient Cable,9ft AHA/IEC, MD6M>LL 5LD, no resistance fit for LL style LeadWires		 105

Compatible Brand	Compatible Models No.	OEM P/N	Evotec Part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
MINDRAY	PM-7000,6201,2004年以后的PM8000,PM9000. with 6-pin input connector	0010-30-12377	EEJ021-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6pin>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires	007 	242 
			EEJ021-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6pin>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires	007 	243 
	PM-6800,2004年以后的PM-6000, with 12-pin green input connector	0010-30-42719	EEJ051-3A/I	3Ld Patient Cable,9ft AHA/IEC, straight Round 12P>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires	222 	242 
			EEJ051-5A/I	5Ld Patient Cable,9ft AHA/IEC, straight Round 12P>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires	222 	243 
MINDRAY	PM6000, with 12-pin green input connector	0010-30-12261	EEJ050-3A/I	3Ld Patient Cable,9ft AHA/IEC, Type L Round 12P>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires	030 	242 
			EEJ050-5A/I	5Ld Patient Cable,9ft AHA/IEC, Type L Round 12P>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires	030 	243 
NEC	3M12	47524	EEK053-3I	3Ld Patient Cable,9ft AHA/IEC, Round 16p>NEC 3LD, no resistance fit for 47504 LeadWires	024 	096 
		47525	EEK053-5I	5Ld Patient Cable,9ft AHA/IEC, Round 16p>NEC 5LD, no resistance fit for 47505 LeadWires	024 	099 
Nihon Kohden	OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K BSM-7105K		EEC037-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 8p>Din 3LD, no resistance fit for Din style LeadWires	237 	101 
		JC-005P	EEN037-3I	3Ld Patient Cable,9ft IEC, Round 8p>NK 6J, no resistance fit for BR-019 LeadWires	395 	155 
	BSM-7102A,OMP-7201C, OMP-7203C, TEC-5200A (s/n<20421),		EEC038-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 11p>Din 3LD, no resistance fit for Din style LeadWires	154 	101 
			EEN038-3I	3Ld Patient Cable,9ft IEC, Round 11p>NK 6J no resistance fit for BR-019 LeadWires	154 	155 

Compatible Brand	Compatible Models No.	OEM P/N	Evotec part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
Nihon Kohden	BSM-2301, BSM-2353, BSM-5100, Multi-parameter bedside monitors	JC-906P	EEP055-6I	6Ld Patient Cable,9ft IEC, Rectangle 12p>NK 8J, no resistance fit for BR-903P, BR-906P	173 	174 
	BSM-8301A/8302A/8500A,BSM-8502A, OEC-6102A(s/n>20156), OEC-6105A(s/n>20226),OEC-7102A/8108A,TEC-5200A(s/n>20421), 7100A/7200A/7300A;	L152 JC-762V TEC-	EEP002-5I	5Ld Patient Cable,9ft IEC, Round 6pin>NK 6J no resistance fit for BR-019 LeadWires	005 	155 
PHILIPS	78352C,7834C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A	M1500A	EED023-3A	3Ld Patient Cable,9ft AHA, Round 12p>AA-3LD 1KQ resistance fit for AA style LeadWires	018 	086 
		M1510A	EED023-3I	3Ld Patient Cable,9ft IEC, Round 12p>AA-3LD 1KQ resistance fit for AA style LeadWires	018 	086 
	78352C,7834C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A	M1520A	EED023-5A	5Ld Patient Cable,9ft AHA, Round 12p>AA-5LD 1KQ resistance fit for AA style LeadWires	018 	085 
		M1530A	EED023-5I	5Ld Patient Cable,9ft IEC, Round 12p>AA-5LD 1KQ resistance fit for AA style LeadWires	018 	085 
PHILIPS	Philips Instruments with 8-Pin Connector: 77000AC,77020A/AC,77020AR,77025A,77030A,77035A,77040A,78208A,78210A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A/B,78354A,78532A/B,78534A/B/C,M1722A/B,M1723A/B,M2400A,M2406A,M2475B,M3500B	M1580A	EED035-3A	3Ld Patient Cable,9ft AHA, Round 8p>AA-3LD 1KQ resistance fit for AA style LeadWires	019 	086 
		M1590	EED035-3I	3Ld Patient Cable,9ft IEC, Round 8p>AA-3LD 1KQ resistance fit for AA style LeadWires	019 	086 
	Philips Instruments with 8-Pin Connector: 78208A,78210A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A,78354A,78532A/B,78534A/B/C,M1722A/B,M3500B	M1600A	EED035-5A	5Ld Patient Cable,9ft AHA, Round 8p>AA-5LD 1KQ resistance fit for AA style LeadWires	019 	085 
			EED035-5I	5Ld Patient Cable,9ft IEC, Round 8p>AA-5LD 1KQ resistance fit for AA style LeadWires	019 	085 
PHILIPS	M3000A, M3001A, M1001A/B,M1002A/B, 78352C,78354C	M1669A	EEDG073-3AI	3Ld Patient Cable,9ft IEC/AAMI, Round 12p>3-LD 1KQ resistance fit for Philips LeadWires		
PHILIPS	M3000A, M3001A, M1001A/B,M1002A/B, 78352C,78354C	M1668A	EEDG073-5AI	5Ld Patient Cable,9ft IEC/AAMI, Round 12p>5-LD 1KQ resistance fit for Philips LeadWires		

OEM	OEM MODELS NUMBER	OEM P/N	Evotec part number	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)
PHILIPS	78834C 78850A	M1620A	EED036-3A	3Ld Patient Cable,9ft AHA, Round 5p>AA-3LD 1KQ resistance fit for AA style LeadWires	236 	086 
			EED036-3I	3Ld Patient Cable,9ft IEC, Round 5p>AA-3LD 1KQ resistance fit for AA style LeadWires	236 	086 
S&W (Artema)	8000 System, Diascope1,2 , Neoscope, DMS 700, defibrillator, Diascope-NT, Athena Systems, Cardio-Aid, DS Traveller, DSNT, DS1,DS2, MCX/LS		EEB018-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires	008 	107 
			EEB018-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires	008 	105 
SPACELABS	90496 Ultraview, 90369 Ultraview,90367 Ultraview	IEC: 700- 0008-07 AHA: 700- 0008-06	EEL009-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 17p>TRU-LINK 5Ld, 1KQ resistance fit for TRU-LINK LeadWires	022 	091 
	GE MEDICAL:Corometrics 503,Corometrics 504,Corometrics 508; SIEMENS:806180A Nuclear Camera		EEJ0136-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 5pin>Euro 3LD, no resistance fit for Euro style LeadWires	236 	242 
SPACELABS BURDICK Viggo-Spectramed (Statham) WESTINGHOUSE	ADVANCED MEDICAL RESEARCH (AMR): Model 3R; AQUITRON MEDICAL: 9101,9216 (Neonate), ATL ULTRASOUND: Ultramark 8,Ultramark 6,Mark VIII,Mark III,Mark II,600,400; BARD: 500,500 GILSON; BECKMAN (Offner): Offner Dynograph (all),VSM-10C,RM-10 BIOCOM; INC: All models; BOURNS LIFE SYSTEMS: 118-1/2/3; DEL MAR: 2000 Series,2820,2920,2920A; GE MEDICAL SYSTEMS: Critikon Exerstress; MEDICAL ELECTRONICS: All Models (except A30, A40); HEALTHWATCH (Cambridge): MC-4,Simpli Scribe, Versa Scribe, VS-1/2/3; INTERSPEC (Irex) (J&J): System II Echo; KDC: All Models (except Infant IM); MARCONI MEDICAL: 464B, 51-690, CA101, CA105, Echoview 80C, ECO-6, RC-101/M; PIONEER (American Optical): Life Rite 50 ; PPG (Datamedix, Becton-Dickinson): 105 w/lead select, ELD-320W, ELD-400, HR-4, IPA-1/2/3/4, IPA-BD21, IS-5LB, MS-25L, S-5BL , Statscope; PPG (E for M): EAI Series,ECG-6 Series,ECI Series,ECV Series,EDI Series,EEP-8,ERI Series,ETD-8 PPG (Honeywell): 5100 w/ lead select ; SPACELABS BURDICK: CS615,CS625A,CSS20 ,CSS23 ,CSS27,CSS7200,EK1,EK2,EK3,EK4,EK5/A/AS,EK8SS,M200,Medic 1,Medic 2,PA-100,PA-104,PA-105,1600,2100; VIGGO-SPECTRAMED (Statham): 6703,Ekoline 21, WESTINGHOUSE: All Models					
Fukuda Denshi	DATASCOPE:Expert FUKUDA DENSHI:DS5100E, DS5300W, DS5400, DS7101/L, DS7141		EEY068-5A/I	5Ld and 3Ld Patient Cable,9ft AHA/IEC, Round 12pin> 16 jack non- resistance fit for EY032S3I LeadWires	253 	373 
Lohmeier			EEC078-3A/I	3Ld Patient Cable,9ft IEC/AAMI , 5pin+5Jack plug>Din-3LD non-resistance fit for Din1.5 style Leads		101 
Lohmeier	D-81241, original P/N.: 98.060.120/A	98.060.1 20/A	EEJM079-5A/I	5Ld Multi-Cable,9ft IEC/AAMI , 27pin plug>Euro 5Ld ECG+2*IBP+SpO2+2*Temp (discharge inside)		

LEAD WIRES


Low profile, insert molded snap or pinch connectors offer patient comfort and secure attachment to snap electrodes. The snap terminatic ergonomic, teardrop design reduces tangling at the patient end connector. Grabber and mini clip terminations have dual sided metal con to provide the connection.

Compatible Brand	Productor Guide	OEM P/N	Evotec part number	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination
Bionet Datex Siemens GE-Hellige Kontron Goldway S&W (Artema)	suitable for all LL-Plug System , Yoke type: Lifelinc(LL),,Code: <u>B</u>		EEB024S3A	LL lead wires, <u>24 in. (61 cm), 3-Ld</u>	AHA (AAMI)	Snap
			EEB024C3A		AHA (AAMI)	Grabber (Clip)
			EEB024S3I		IEC	Snap
			EEB024C3I		IEC	Grabber (Clip)
			EEB024S5A	LL lead wires, <u>24 in. (61 cm), 5-Ld</u>	AHA (AAMI)	Snap
			EEB024C5A		AHA (AAMI)	Grabber (Clip)
			EEB024S5I		IEC	Snap
			EEB024C5I		IEC	Grabber (Clip)
			EEB040S3A	LL lead wires, <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Snap
			EEB040C3A		AHA (AAMI)	Grabber (Clip)
			EEB040S3I		IEC	Snap
			EEB040C3I		IEC	Grabber (Clip)
			EEB040S5A	LL lead wires, <u>40 in. (102 cm), 5-Ld</u>	AHA (AAMI)	Snap
			EEB040C5A		AHA (AAMI)	Grabber (Clip)
			EEB040S5I		IEC	Snap
			EEB040C5I		IEC	Grabber (Clip)
Biosys Creative Datascope Nihon Kohden AIR SHIELDS (NARCO, Vickers) ATL ULTRASOUND BIOCHEM INTERNATIONAL COLIN Medical CRITICARE IVY	suitable for all Din-Plug System , Yoke type: Din style,Code: <u>C</u>		EEB018S3A	Din lead wires, <u>18 in. (45 cm), 3-Ld</u>	AHA (AAMI)	Snap
			EEB018C3A		AHA (AAMI)	Grabber (Clip)
			EEC018S3I		IEC	Snap
			EEC018C3I		IEC	Grabber (Clip)
			EEC018S5A	Din lead wires, <u>18 in. (45 cm), 5-Ld</u>	AHA (AAMI)	Snap
			EEC018C5A		AHA (AAMI)	Grabber (Clip)
			EEC018S5I		IEC	Snap
			EEC018C5I		IEC	Grabber (Clip)
			EEC024S3A	Din lead wires, <u>24 in. (61 cm), 3-Ld</u>	AHA (AAMI)	Snap
			EEC024C3A		AHA (AAMI)	Grabber (Clip)
			EEC024S3I		IEC	Snap
			EEC024C3I		IEC	Grabber (Clip)
			EEC024S5A	Din lead wires, <u>24 in. (61 cm), 5-Ld</u>	AHA (AAMI)	Snap
			EEC024C5A		AHA (AAMI)	Grabber (Clip)
			EEC024S5I		IEC	Snap
			EEC024C5I		IEC	Grabber (Clip)
			EEC040S3A	Din lead wires, <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Snap
			EEC040C3A		AHA (AAMI)	Grabber (Clip)
			EEC040S3I		IEC	Snap
			EEC040C3I		IEC	Grabber (Clip)
			EEC040S5A	Din lead wires, <u>40 in. (102 cm), 5-Ld</u>	AHA (AAMI)	Snap
			EEC040C5A		AHA (AAMI)	Grabber (Clip)
			EEC040S5I		IEC	Snap
			EEC040C5I		IEC	Grabber (Clip)

Compatible Brand	Product Guide	OEM P/N	Evotec part number	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination
PHILIPS M3000A, M3001A, M1001A/B, M1002A/B, 78352C, 78354C		M1673A	EEDG040S3A	Multi-Link Ldwr, <u>3-Ld</u> Set - Group, <u>40 in (102cm)</u>	AHA (AAMI)	Snap
		M1674A	EEDG040S3I		IEC	Snap
		M1671A	EEDG040C3A		AHA (AAMI)	Grabber (Clip)
		M1672A	EEDG040C3I		IEC	Grabber (Clip)
		M1644A	EEDG040S5A	Multi-Link Ldwr, <u>5-Ld</u> Set - Group, LA, RA, V=40in(102cm); LL, RL=61in(155cm);	AHA (AAMI)	Snap
		M1645A	EEDG040S5I		IEC	Snap
		M1968A	EEDG040C5A		AHA (AAMI)	Grabber (Clip)
		M1971A	EEDG040C5I		IEC	Grabber (Clip)
PHILIPS  139	 ED040C5A	M1605A	EED040S3A	AA lead wires, <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Snap
		M1603A	EED040C3A		AHA (AAMI)	Grabber (Clip)
		M1615A	EED040S3I		IEC	Snap
		M1613A	EED040C3I		IEC	Grabber (Clip)
		M1625A	EED040S5A	AA lead wires, <u>5-Ld</u> LL, RL, V= <u>40in. (102cm)</u> , <u>RA, RL=51in. (130 cm)</u>	AHA (AAMI)	Snap
		M1623A	EED040C5A		AHA (AAMI)	Grabber (Clip)
		M1635A	EED040S5I		IEC	Snap
		M1633A	EED040C5I		IEC	Grabber (Clip)
		M1601A	EED040C3A-OR	AA lead wires, <u>OR</u> <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Grabber (Clip)
		M1611A	EED040C3I-OR		IEC	Grabber (Clip)
		M1621A	EED040C5A-OR	AA lead wires, <u>5-Ld, OR</u> LL, RL, V=40in. (102cm), RA, RL=51in. (130 cm)	AHA (AAMI)	Grabber (Clip)
		M1631A	EED040C5I-OR		IEC	Grabber (Clip)
GE--Marquette (Set -Group type Lead wires)  185	 EE029C3A	E9008LF	EEE029S3A-G	Multi-Link Ldwr Set - Group,, <u>29 in. (76 cm), 3-Ld</u>	AHA (AAMI)	Snap
		E9003CL	EEE029C3A-G		AHA (AAMI)	Grabber (Clip)
		E9008LH	EEE029S3I-G		IEC	Snap
		E9003CN	EEE029C3I-G		IEC	Grabber (Clip)
		E9008LG	EEE051S3A-G	Multi-Link Ldwr Set - Group,, <u>51 in. (130 cm), 3-Ld</u>	AHA (AAMI)	Snap
		E9003CM	EEE051C3A-G		AHA (AAMI)	Grabber (Clip)
		E9008LJ	EEE051S3I-G		IEC	Snap
		E9003CP	EEE051C3I-G		IEC	Grabber (Clip)
GE--Marquette (Set -Group type Lead wires) 	 EE029S5A	E9008L	EEE029S5A-G	Multi-Link Ldwr Set - Group,, <u>29 in. (76 cm), 5-Ld</u>	AHA (AAMI)	Snap
		E9003CE	EEE029S5A-G		AHA (AAMI)	Grabber (Clip)
		E9008LB	EEE029S5I-G		IEC	Snap
		E9003CG	EEE029C5I-G		IEC	Grabber (Clip)
		E9008LA	EEE051S5A-G	Multi-Link Ldwr Set - Group,, <u>51 in. (130 cm), 5-Ld</u>	AHA (AAMI)	Snap
		E9003CF	EEE051C5A-G		AHA (AAMI)	Grabber (Clip)
		E9008LC	EEE051S5I-G		IEC	Snap
			EEE051C5I-G		IEC	Grabber (Clip)
		E9008LD	EEE051S5A-GM	Multi-Link Ldwr SetGroup, <u>5-Ld, Mixed Length,</u> RL,LL,V=51in.(130cm) ; RA,LA=29 in.(74cm)	AHA (AAMI)	Snap
		E9003CJ	EEE051C5A-GM		AHA (AAMI)	Grabber (Clip)
		E9008LE	EEE051S5I-GM	Multi-Link Ldwr SetGroup, <u>5-Ld, Mixed Length,</u> N,F,C=51in. (130cm) , R,L=29 in.(74cm)	IEC	Snap
			EEE051C5I-GM		IEC	Grabber (Clip)
GE-Medical Marquette (Set -Assembly type Lead wires)  140	 EE051C5A-M	E9008KB	EEE029S5A	Multi-Link Ldwr Set <u>29 in. (76 cm), 5-Ld</u>	AHA (AAMI)	Snap
		E9002ZW	EEE029C5A		AHA (AAMI)	Grabber (Clip)
		E9008KD	EEE029S5I		IEC	Snap
		E9002ZZ	EEE029C5I		IEC	Grabber (Clip)
		E9008KC	EEE051S5A	Multi-Link Ldwr Set, <u>51 in. (130 cm), 5-Ld</u>	AHA (AAMI)	Snap
		E9002ZY	EEE051C5A		AHA (AAMI)	Grabber (Clip)
		E9008KE	EEE051S5I		IEC	Snap
		E9003A	EEE051C5I		IEC	Grabber (Clip)
		E9008KF	EEE051S5A-M	Multi-Link Ldwr Set, <u>5-Ld, Mixed Length,</u> RL,LL,V=51in.(130cm); RA,LA=29 in.(74cm)	AHA (AAMI)	Snap
		E9003AA	EEE051C5A-M		AHA (AAMI)	Grabber (Clip)
			EEE051S5I-M	Multi-Link Ldwr Set, <u>5-Ld, Mixed Length,</u> N,F,C=51in.(130cm) ; R,L=29 in.(74cm)	IEC	Snap
		E9003AB	EEE051C5I-M		IEC	Grabber (Clip)

Compatible Brand	Product Guide	OEM P/N	Evotec part number	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination
Datex PHILIPS PIONEER PPG SPACELABS BURDICK TEKTRONIX (Vitatek) WELCH ALLYN ZOLL  138	suitable for all Conmed-Plug System , Yoke type: <u>G</u>  EG040S5A		EEG024S3A		AHA (AAMI)	Snap
			EEG024C3A	Conmed lead wires,	AHA (AAMI)	Grabber (Clip)
			EEG024S3I	<i>24 in. (61 cm), 3-Ld</i>	IEC	Snap
			EEG024C3I		IEC	Grabber (Clip)
			EEG024S5A		AHA (AAMI)	Snap
			EEG024C5A	Conmed lead wires,	AHA (AAMI)	Grabber (Clip)
			EEG024S5I	<i>24 in. (61 cm), 5-Ld</i>	IEC	Snap
			EEG024C5I		IEC	Grabber (Clip)
			EEG040S3A		AHA (AAMI)	Snap
			EEG040C3A	Conmed lead wires,	AHA (AAMI)	Grabber (Clip)
			EEG040S3I	<i>40 in. (102 cm), 3-Ld</i>	IEC	Snap
			EEG040C3I		IEC	Grabber (Clip)
			EEG040S5A		AHA (AAMI)	Snap
			EEG040C5A	Conmed lead wires,	AHA (AAMI)	Grabber (Clip)
			EEG040S5I	<i>40 in. (102 cm), 5-Ld</i>	IEC	Snap
			EEG040C5I		IEC	Grabber (Clip)
Datex-Ohmeda GE-medical Pro1000  137	suitable for all ASP-Plug System , Yoke type: <u>H</u>  EH024C3I		EEH024S3A		AHA (AAMI)	Snap
			EEH024C3A	ASP lead wires,	AHA (AAMI)	Grabber (Clip)
			EEH024S3I	<i>24 in. (61 cm), 3-Ld</i>	IEC	Snap
			EEH024C3I		IEC	Grabber (Clip)
			EEH024S5A		AHA (AAMI)	Snap
			EEH024C5A	ASP lead wires,	AHA (AAMI)	Grabber (Clip)
			EEH024S5I	<i>24 in. (61 cm), 5-Ld</i>	IEC	Snap
			EEH024C5I		IEC	Grabber (Clip)
			EEH040S3A		AHA (AAMI)	Snap
			EEH040C3A	ASP lead wires,	AHA (AAMI)	Grabber (Clip)
			EEH040S3I	<i>40 in. (102 cm), 3-Ld</i>	IEC	Snap
			EEH040C3I		IEC	Grabber (Clip)
			EEH040S5A		AHA (AAMI)	Snap
			EEH040C5A	ASP lead wires,	AHA (AAMI)	Grabber (Clip)
			EEH040S5I	<i>40 in. (102 cm), 5-Ld</i>	IEC	Snap
			EEH040C5I		IEC	Grabber (Clip)
type LM Lead wires  494			EEZ040S3A		AHA (AAMI)	Snap
			EEZ040C3A	LM lead wires,	AHA (AAMI)	Grabber (Clip)
			EEZ040S3I	<i>40 in. (102 cm), 3-Ld</i>	IEC	Snap
			EEZ040C3I		IEC	Grabber (Clip)
			EEZ040S5A		AHA (AAMI)	Snap
			EEZ040C5A	LM lead wires,	AHA (AAMI)	Grabber (Clip)
			EEZ040S5I	<i>40 in. (102 cm), 5-Ld</i>	IEC	Snap
			EEZ040C5I		IEC	Grabber (Clip)
Mindray Siemens SPACELABS BURDICK Viggo-Spectramed (Statham) WESTINGHOUSE  084	suitable for all Euro-Plug System , Yoke type: <u>J</u>  EJ040S3I  EJ040S5I		EEJ024S3A		AHA (AAMI)	Snap
			EEJ024C3A	Euro style lead wires,	AHA (AAMI)	Grabber (Clip)
			EEJ024S3I	<i>24 in. (61 cm), 3-Ld</i>	IEC	Snap
			EEJ024C3I		IEC	Grabber (Clip)
			EEJ024S5A		AHA (AAMI)	Snap
			EEJ024C5A	Euro style lead wires,	AHA (AAMI)	Grabber (Clip)
			EEJ024S5I	<i>24 in. (61 cm), 5-Ld</i>	IEC	Snap
			EEJ024C5I		IEC	Grabber (Clip)
			EEJ040S3A		AHA (AAMI)	Snap
			EEJ040C3A	Euro style lead wires,	AHA (AAMI)	Grabber (Clip)
			EEJ040S3I	<i>40 in. (102 cm), 3-Ld</i>	IEC	Snap
			EEJ040C3I		IEC	Grabber (Clip)
			EEJ040S5A		AHA (AAMI)	Snap
			EEJ040C5A	Euro style lead wires,	AHA (AAMI)	Grabber (Clip)
			EEJ040S5I	<i>40 in. (102 cm), 5-Ld</i>	IEC	Snap
			EEJ040S3A		IEC	Grabber (Clip)

Compatible Brand	Productor Guide	OEM P/N	Evotec part number	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination
NEC suitable for all NEC-Plug System , Yoke type: NEC style, Code: K	 097		EEK029S3I	47504 style lead wires, <u>29 in. (76 cm), 3-Ld</u>	IEC	Snap
		47504	EEK029C3I		IEC	Grabber (Clip)
	 100		EEK029S5I	47505 style lead wires, <u>29 in. (76 cm), 5-Ld</u>	IEC	Snap
		47505	EEK029C5I		IEC	Grabber (Clip)
SPACELABS  166	suitable for all TRU-LINK Plug System , Yoke type: TRU-LINK style, Code: <u>L</u>	700-0007-00	EEL024S3A	TUR-LINK style lead wires, <u>24 in. (61 cm), 3-Ld</u>	AHA (AAMI)	Snap
		700-0008-00	EEL024C3A		AHA (AAMI)	Grabber (Clip)
		700-0007-01	EEL024S3I		IEC	Snap
		700-0008-01	EEL024C3I		IEC	Grabber (Clip)
		700-0007-04	EEL024S4A	TUR-LINK style lead wires, <u>24 in. (61 cm), 4-Ld</u>	AHA (AAMI)	Snap
		700-0008-04	EEL024C4A		AHA (AAMI)	Grabber (Clip)
		700-0007-05	EEL024S4I		IEC	Snap
		700-0008-05	EEL024C4I		IEC	Grabber (Clip)
		700-0007-08	EEL024S5A	TUR-LINK style lead wires, <u>24 in. (61 cm), 5-Ld</u>	AHA (AAMI)	Snap
		700-0008-08	EEL024C5A		AHA (AAMI)	Grabber (Clip)
		700-0007-09	EEL024S5I		IEC	Snap
		700-0008-09	EEL024C5I		IEC	Grabber (Clip)
		700-0007-12	EEL024S6A	TUR-LINK style lead wires, <u>24 in. (61 cm), 6-Ld</u>	AHA (AAMI)	Snap
		700-0008-12	EEL024C6A		AHA (AAMI)	Grabber (Clip)
		700-0007-13	EEL024S6I		IEC	Snap
		700-0008-13	EEL024C6I		IEC	Grabber (Clip)
		700-0007-16	EEL024S10A	TUR-LINK style lead wires, <u>24 in. (61 cm), 10-Ld</u>	AHA (AAMI)	Snap
		700-0008-16	EEL024C10A		AHA (AAMI)	Grabber (Clip)
		700-0007-17	EEL024S10I		IEC	Snap
		700-0008-17	EEL024S10I		IEC	Grabber (Clip)
		700-0007-02	EEL040S3A	TUR-LINK style lead wires, <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Snap
		700-0008-02	EEL040C3A		AHA (AAMI)	Grabber (Clip)
		700-0007-03	EEL040S3I		IEC	Snap
		700-0008-03	EEL040C3I		IEC	Grabber (Clip)
		700-0007-06	EEL040S4A	TUR-LINK style lead wires, <u>40 in. (102 cm), 4-Ld</u>	AHA (AAMI)	Snap
		700-0008-06	EEL040C4A		AHA (AAMI)	Grabber (Clip)
		700-0007-07	EEL040S4I		IEC	Snap
		700-0008-07	EEL040C4I		IEC	Grabber (Clip)
		700-0007-10	EEL040S5A	TUR-LINK style lead wires, <u>40 in. (102 cm), 5-Ld</u>	AHA (AAMI)	Snap
		700-0008-10	EEL040C5A		AHA (AAMI)	Grabber (Clip)
		700-0007-11	EEL040S5I		IEC	Snap
		700-0008-11	EEL040C5I		IEC	Grabber (Clip)
		700-0007-14	EEL040S6A	TUR-LINK style lead wires, <u>40 in. (102 cm), 6-Ld</u>	AHA (AAMI)	Snap
		700-0008-14	EEL040C6A		AHA (AAMI)	Grabber (Clip)
		700-0007-15	EEL040S6I		IEC	Snap
		700-0008-15	EEL040C6I		IEC	Grabber (Clip)
		700-0007-18	EEL040S10A	TUR-LINK style lead wires, <u>40 in. (102 cm), 10-Ld</u>	AHA (AAMI)	Snap
		700-0008-18	EEL040C10A		AHA (AAMI)	Grabber (Clip)
		700-0007-19	EEL040S10I		IEC	Snap
		700-0008-19	EEL040C10I		IEC	Grabber (Clip)

Compatible Brand	Productor Guide	OEM P/N	Evotec part number	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination
Colin suitable for all Colin-Plug System , Yoke type: Colin style, Code: <u>M</u>	 009	CL-3601	EEM040S3I	COLIN style lead wires, <u>40in. (102 cm), 3-Ld</u>	IEC	Snap
			EEM040C3I		IEC	Grabber (Clip)
	 010		EEM040S5I	COLIN style lead wires, <u>40in. (102 cm), 5-Ld</u>	IEC	Snap
			EEM040C5I		IEC	Grabber (Clip)
Nihon Konden suitable for all NK-5LD Plug System , Yoke type: NK-5LD style, Code: <u>N</u>	 126	BR-019	EEN032S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap
		BR-019P	EEN032C3I		IEC	Grabber (Clip)
	 125	BR-021P	EEN032S5I	NK style lead wires, <u>32 in. 80 cm, 5-Ld</u>	IEC	Snap
			EEN032C5I		IEC	Grabber (Clip)
Nihon Konden suitable for all NK-6LD Plug System , Yoke type: NK-6LD style, Code: <u>P</u>	 171/433	BR-913P	EEP035S3I	NK style lead wires, <u>35 in. 90 cm, 3-Ld</u>	IEC	Snap
		BR-903P	EEP035C3I		IEC	Grabber (Clip)
			EEP035S6I	NK style lead wires, <u>35 in. 90 cm, 6-Ld</u>	IEC	Snap
		BR-906P	EEP035C6I		IEC	Grabber (Clip)
Nihon Konden suitable for all NK-5LD Plug System , Yoke type: NK-5LD style, Code: <u>R</u>	 124		EHA004S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap
		BR-546S	EHA004C3I		IEC	Grabber (Clip)
	 123		EEV032S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap
		BR-002P	EEV032C3I		IEC	Grabber (Clip)
	 122		EEV032S5I	NK style lead wires, <u>32 in. 80 cm, 5-Ld</u>	IEC	Snap
		BR-004P	EEV032C5I		IEC	Grabber (Clip)
FUKUDA DENSHI	 372		EEY032S3I	one-piece lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap
			EEY032S3A		AHA (AAMI)	Snap



EED040C3A



EED040S3



EEE029S3A-G



EHA002S



EEC040S3A