

GD 300

Electric surgical unit

Guoteng

**One Of The Professional Chinese
MANUFACTURERS
Special In
Medical Equipment**



Features:

The wide scope permanent power, the surgery locates precisely.


The return route pole plate monitor system, guarantees the patient safety

Seven kind of working model, the surgery is quickly easy.

(pure cutting, blend 1, blend 2, coagulation, spay, bipolar cut, bipolar coagulation)

Many kinds of appendix disposition, best performance-to-price ratio

All the connector compatible with international standard.



Guoteng Science & Technology Development Co., Ltd. Zhuhai

Add.: 6/F., Three-Dimension Science & Technology Building, No. 184 Bailian Road, Jida, Zhuhai, Guangdong, 519015 China

Tel: +86-756-3321860 Fax: +86-756-6181801 E-mail: market@guoteng.gd.cn

<http://www.guoteng.gd.cn>

Specification:

		Work Mode	Current Intensity	Rated Power	Wave	Vp-p
Monopole	Incision	Pure	400KHz	0—300W (400 Ω)	Sine Wave	2000V
		Blend I	400KHz	0—250W (400 Ω)	Modulate Sine Wave	3000V
		Blend II	400KHz	0—200W (400 Ω)	Modulate Sine Wave	3700V
	Concreting	Standard	400KHz	0—100W (400 Ω)	Modulate Sine Wave	4700V
		Spray	400KHz	0—80W (400 Ω)	Modulate Sine Wave	6000V
Bipolar	Coagulation	Cut	400KHz	0—100W (100 Ω)	Sine Wave	1000V
		Coag	400KHz	0—70W (100 Ω)	Modulate Sine Wave	700V

Application

Ordinary surgery, cardiology, urology, gynaecology, chest chirurgery, oncology

ACCESSORIES FULL FIT FOR THE INCISION AND MELTING OF THE PARENCHYMA**Other Guoteng Products:**

**GR 120 Radio
electric wave bistoury**

Application

Gynaecologic, the five sense organs and neurosurgery, ENT, and cosmetology.

Guoteng Science & Technology Development Co., Ltd. Zhuhai

Add.: 6/F., Three-Dimension Science & Technology Building, No. 184 Bailian Road, Jida, Zhuhai, Guangdong, 519015 China

Tel: +86-756-3321860 Fax: +86-756-6181801 E-mail: market@guoteng.gd.cn

<http://www.guoteng.gd.cn>