

	Monitor / Recorder	Manufacturer Chart #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Advanced Medical	IM 88	1142752-012A/ 31-0012	7778108012	CTG	Pack	120 x 100 x 150	Red	5
	IM76/66	1142751-012E	31-0410	CTG	Pack	120 x 100 x 150	Red	5
	IM77	31-0427	31-0427	CTG	Pack	150 x 100 x 150	Brown	5
Ames	CLINITEK 50			LAB	Roll	57 x 30	White	5
	CLINITEK 100	5773	5773	LAB	Roll	57 x 30	White	5
Analogic	U-S-A/INT-3	82275	82275	CTG	Pack	152 x 90 x 160	Red	5
		CREATIVE PC800	CREATIVE PC800	CTG	Pack	110x100x150	Red	5
Antinatal Fetal	Fetaltrack 310	FM1D122-101	FM1D122-101	CTG	Pack	114 x 100 x 150	Blue	1
Aspel	B5			ECG	Roll	58 x 25	Orange	5
	B56			ECG	Roll	112 x 25	Red	5
Biocare	300			ECG	Roll	63 x 30	Orange	10
	1200			ECG	Roll	210 x 30	Orange	5
		Biocare 801		CTG	Pack	112x100x160	Orange	5
		101G	101G	ECG	Roll	50 x 30	Green	10
Bionet	Cardiocare	0700-01	0700-01	CTG	Roll	215 x 46OD	Green	1
		FC-700	FC-700	ECG	Roll	214 x 30	Orange	1
		FC-1400	FC-1400	CTG	Roll	214 x 30	Red	1
		BFM900L	BFM900L	CTG	Roll	151 x 30	Green	1
Bionics				CTG	Roll	111.5x48OD	Orange	5
Biosys	BFM800			CTG	Pack	112 x 120 x 250	Blue	5
	IFM500	94-001	94-001	CTG	Pack	120 x 120 x 250	Blue	5
Bistos	AFM-0500	FS130-120-30R-01	FS130-120-30R-01	CTG	Pack	130 x 120 x 250	Red	1
	MF400/BT300	FS130-120-30R	FS130-120-30R	CTG	Pack	130 x 120 x 250	Red	1
	BT-350	FS151-90-80R-01	FS151-90-80R-01	CTG	Pack	152 x 90 x 150	Red	5
	BTL-08	BTL-08	BTL08	SPIRO	Roll	111.5 x 25	Red	5
Burdick	EK 10/ELITE II	007957/58	007957/58	ECG	Roll	50 x 30	Red	10
	ECLIPSE 4/ ECLIPSE 400	007970/007971	007970/007971	ECG	Pack	106 x 140 x 200	Red	1
	ECLIPSE 8	007986	007986	ECG	Pack	210 x 300 x 200	Orange	1
	E350	007868	007868	ECG	Pack	215 x 280 x 215	Red	1
	ECLIPSE 850	007980	007980	ECG	Pack	210 x 150 x 300	Orange	1
Cadence	MFM-2	MSI-01921	MSI-01921	CTG	Pack	112 x 90 x 150	Red	5
	MFM2	M50-78019	M50-78019	CTG	Pack	112x90x150	Red	5
Cardioline	DELTA 1	8743550	8743550	ECG	Roll	50 x 30	Red	10
	DELTA 1 Plus/i	8743560	8743560	ECG	Roll	60 x 30	Red	10
	DELTA 1 Plus/i	8745558	8745558	ECG	Pack	60 x 100 x 300	Red	5
	DELTA 3	8745555	8745555	ECG	Pack	134 x100 x 300	Orange	5
	DELTA 3 Plus/6i	8745554	8745554	ECG	Pack	112 x 100 x 300	Orange	5
	DELTA 30	8743555	8743555	ECG	Roll	183 x 30	Orange	1
	DELTA 60	8743555	8743555	ECG	Pack	183 x 270 x 200	Red	1
	DELTA 60 Plus	8743556	8743556	ECG	Roll	210 x 30	Orange	1
	DELTA 612	8745556	8745556	ECG	Pack	210 x 280 x 215	Red	1

	Monitor / Recorder	Manufacturer #	Chart	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Cardiorapid	K111 / DOC101 K312	60100040		60100040	ECG	Roll	50 x 22	Red	5
					ECG	Pack	50 x 70 x 200	Red	5
Carewell		1103G Carewell 6 channel Carewell 12 channel		1103G Carewell 6 channel Carewell 12 channel	ECG	Roll	63 x 30	Orange	10
					ECG	Roll	112x45OD	Red	5
Ciba Corning	600 SERIES EXPRESS 560 200 SERIES M800 M102	478638 233720 473550 673709000 478368		478638 233720 473550 673709000 478368	ECG	Roll	216x45OD	Red	5
					LAB	Roll	44 x 51OD	Bond	2
Cobas Mira	HUMALYSER 2000 PRESS-MATE 8800P PULSEMATE BX50	1021125		1021125	LAB	Roll	110 x 50OD	White	5
					LAB	Roll	62 x 73OD	White	5
Collin					LAB	Roll	110 x 47OD	White	5
					LAB	Roll	60 x 15	Silver	5
Collins	SPIROMETER				LAB	Roll	140 x 36	White	4
					LAB	Roll	57 x 30	White	5
Corometrics	115-150 170 151 171 145 145	4305AAO 4305BAO 4305CAO 4305DAO 4483AAO 4483BAO		4305AAO 4305BAO 4305CAO 4305DAO 4483AAO 4483BAO	BP	Roll	57 x 25	White	5
					BP	Roll	57 x 47OD	White	5
Cosmed	COSMED PONY NIBP 1100 NIBP 1200				SPIRO	Roll	58 x 25	White	5
					CTG	Roll	152 x 90 x 150	Red	5
Critikare	WAVEFORM PRO 100/300/400 SXP DINAMAP IR2000 COMPACT IS	8730 1846 8214-8955 8255 8956		8730 1846 8214-8955 8255 8956	CTG	Pack	152 x 90 x 150	Red	5
					CTG	Pack	152 x 90 x 150	Red	5
CMS800G	310B 310B	CMS800G 49030270 49020130		49030270 49020130	CTG	Roll	140 x 42	Red	1
					CTG	Roll	140 x 42	Red	1
Dantec	URODYN 100 3000/3000A/PASSPORT EL BALLON PUMP/SYSTEMS 90/95	9021B3111 0683-00-0422-01/02 0063-00-0241-02		9021B3111 0683-00-0422-01/02 0063-00-0241-02	SPIRO	Roll	57 x 10	Bond	5
					BP	Roll	58 x 30	White	5
Datascop	Cardicap 11	74202 73021		74202 73021	BP	Roll	50 x 30	White	10
					BP	Roll	50 x 30	White	10
Datex	SE1200 EDAN F9-G EDAN F9-G	M25R-7511 M25R-7512		M25R-7511 M25R-7512	MED	Roll	50 x 25	White	5
					MED	Roll	50 x 25	White	5
Edan Instruments		MSI-19901 Edan SE3 MSI-32462/ 9401-021		MSI-19901 Edan SE3 MSI-32462/ 9401-021	CTG	Pack	152x90x150	Pink	5
					CTG	Pack	151x90x150	Green	5

	Monitor / Recorder	Manufacturer Chart #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Electronica Trentina	DOT-CARD/AUTORULER 12/14/DOC101	66010001/A/B	66010001/A/B	ECG	Roll	50 x 22	Red	5
	DOT-CARD	66010002/B	66010002/B	ECG	Roll	134 x 25	Blue	5
	EXCELL 106-103	66010008	66010008	ECG	Roll	120 x 25	Blue	5
	START 100	66010026	66010026	ECG	Roll	130 x 25	Blue	5
	AR1200	66010033	66010033	ECG	Roll	120 x 20	Blue	5
	AR600	66010037	66010037	ECG	Roll	60 x 18	Blue	5
	Start 200 / Daedalus 346	66010029	66010029	ECG	Roll	210 x 20	Blue	5
	AR2100ADV	66010027	66010027	ECG	Pack	208x150x25.5	Blue	1
	Daedalus 346	66010022	66010022	ECG	Pack	210 x 300 x 200	White	1
	HX20			SPIRO	Roll	57 x 10	Bond	5
Epson								
	PERSONAL C1	99220086138	99220086138	ECG	Roll	50 x 30	Green	10
	P80	9690029000	9690029000	ECG	Pack	90 x 70 x 400	Red	5
	C120	9690011000	9690011000	ECG	Pack	126 x 150 x 170	Red	1
	PERSONAL C3	99220098229	99220098229	ECG	Roll	142 x 30	Yellow	5
	ESB5 / C210	9690011001	9690011001	ECG	Pack	210 x 150 x 170	Red	1
	P80	9690029020	9690029020	ECG	Roll	90 x 28	Red	5
	FORMULA	99330017138	99330017138	ECG	Pack	210 x 300 x 200	Red	1
	ACTAVIS	9690027000	9690027000	ECG	Pack	210 x 150 x 200	White	1
	6400/6500			AUD	Roll	57 x 45	White	5
Fonix	6400CX			AUD	Roll	59 x 30	White	5
	FP30/40			AUD	Roll	59 x 30	White	5
Fukuda-Denshi	Cardiosunny 300	OP12KE	OP12KE	ECG	Roll	50 x 30	Orange	10
	Cardiosunny 503	OP30KE	OP30KE	ECG	Roll	130 x 30	Orange	5
	FCP15-G/FCP11	OP110TE	OP110TE	ECG	Roll	63 x 30	Orange	10
	FCP231G/FX7404	OP69TE	OP69TE	ECG	Roll	210 x 30	Orange	5
	FX102	OP115TE	OP115TE	ECG	Pack	57 x 75 x 400	Orange	10
	FCP15\FCP11\FX2111	OP119TE	OP119TE	ECG	Roll	63 x 30	Orange	10
	FX7202	OP222TE	OP222TE	ECG	Pack	110 x 140 x 143	Orange	1
	FCP220G-FX2201X	OP122TE	OP122TE	ECG	Pack	63 x 75 x 400	Orange	10
	FD31-FD35\35\36	OP31KE	OP31KE	ECG	Roll	145 x 30	Orange	5
	FX3101	OP358TE	OP358TE	ECG	Roll	145 x 30	Orange	5
Fukuda ME	Cardiosunny 500/501/501D	050R30-F01	050R30-F01	ECG	Roll	50 x 30	Orange	10
	Cardisunny 501A	050F30-H02	050F30-H02	ECG	Pack	50 x 70 x 400	Orange	10
	ALPHA 1000/FX7101	063R30	063R30	ECG	Roll	63 x 30	Orange	10
	Cardisunny 501AX	063F30-H03	063F30-H03	ECG	Pack	63 x 100 x 300	Orange	10
	Cardisunny 503A/T130	130R30-F01	130R30-F01	ECG	Roll	130 x 30	Orange	5
	503SX-503S	130F50-H02	130F50-H02	ECG	Pack	130 x 75 x 666	Orange	1
Fukuda-Sangyo	SPIROMETER ST95/ST98			SPIRO	Roll	57 x 25	White	5
	ST 250			SPIRO	Roll	110 x 25	White	5
	SP 200			SPIRO	Roll	58 x 25	White	5
	SP3000-SP5000			SPIRO	Roll	110 x 25	White	5

	Monitor / Recorder	Manufacturer #	Chart #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Grason Stadler	GSI 28				AUD	Roll	57 x 65OD	Black	5
	GSI 33/38/MICROTYMP 2				AUD	Roll	110 x 30	White	5
	GSI 67/68/37				AUD	Roll	57 x 25	White	5
	EK33 / 34 / 43	226 061 03		226 061 03	ECG	Roll	130 x 40	Orange	1
	EK41/51	226 130 01		226 130 01	ECG	Roll	55 x 36	Orange	5
Hellige	DEFIB/ CARDIOSERVE	226 130 02		226 130 02	ECG	Roll	55 x 36	White	5
	EK431	226 133 01		226 133 01	ECG	Pack	120 x 150 x 500	Orange	2
	EK512	226 145 02		226 145 02	ECG	Pack	207 x 135 x 370	White	1
	EK53 / 56	226 146 01		226 146 01	ECG	Pack	130 x 135 x 370	Orange	2
	SMR210	226 150 01		226 150 01	ECG	Pack	50 x 100 x 500	White	10
Hellige	SCP900	226 151 01		226 151 01	ECG	Roll	50 x 25	White	10
	MIDAS	226 165 01		226 165 01	ECG	Pack	207 x 280 x 500	White	1
	CARDIOSMART	226 166 01		226 166 01	ECG	Pack	210 x 300 x 150	Orange	1
	MAC 1200	226 166 05		226 166 05	ECG	Pack	210 x 295 x 150	Orange	1
	MICROSMART/MAC500	226 167 02		226 167 02	ECG	Pack	90 x 90 x 388	Orange	1
Hewlett Packard	CONTRAST	226 168 02		226 168 02	ECG	Roll	90 x 35	Orange	1
	43100A\43200A	40453A		40453A	ECG	Roll	50 x 30	Red	10
	Codemaster	40457C/D		40457C/D	ECG	Roll	50 x 30	Grey	10
	M1125A	40477A		40477A	ECG	Roll	50 x 30	White	10
	HP1350A SERIES 50/1351A/1353/2922A	M1911A		M1911A	CTG	Pack	151 x 100 x 150	Green	5
	8025.8030\2.8040.5040\1	9270-0484		9270-0484	CTG	Pack	151 x 100 x 150	Green	5
	Page Writer	4700AH-K21		4700AH-K21	ECG	Page	297 x 210	Red	1000
	Page Writer	4700AH-K22		4700AH-K22	ECG	Page	297 x 210	Black	1000
	M1702A / PW 100 / 200	M1709A		M1709A	ECG	Pack	210 x 300 x 200	Red	8
	M1702A	M1710A		M1710A	ECG	Pack	210 x 300 x 200	Red	8
Honeywell	Pagewriter	M2481A		M2481A	ECG	Pack	216 x 280 x 200	Red	1
	M1702A / PW 100 / 200	M2483A		M2483A	ECG	Pack	210 x 300 x 200	Red	8
	TRIM 1	M3708A		M3708A	ECG	Pack	210 X 300 X 100	Red	1
	M1117A	40469A/B		40469A/B	ECG	Pack	120 x 150 x 300	White	1
	4700 Series	01099415		01099415	ECG	Pen		Black	1
Hokanson	CR100	747 00098		747 00098	ECG	Roll	50 x 20	Orange	5
	Humphrey	TCPI10		TCPI10	MED	Roll	50 x 50OD	White	10
					OPT	Roll	214 x 60OD	White	1
	TECH 15 / Baby Dop	BD 3000/4000		BD 3000/4000	CTG	Pack	150 x 150 x 200	White	1
	HNE	BD4000		BD4000	CTG	Pack	210 x 150 x 200	White	1
ICS Medical	Audiology				AUD	Pack	125 x 150 x 274	White	1
	IFM500 (Refer to Biosys)	94-001		94-001					
	Heart Mirror	HM-3/D		HM-3/D	ECG	Roll	57 x 25	Orange	5
	Heart Screen	HS-3/D		HS-3/D	ECG	Roll	80 x 47OD	Orange	5
		JPD-300A		JPD-300A	CTG	Pack	112 x 100 x 150	Red/Green	5
Jaeger	FLOW MATE/FLOWSCREEN	891-109		891-109	SPIRO	Roll	110 x 30	White	5

	Monitor / Recorder	Manufacturer #	Chart #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Kenz	103/106	R50X30R-C		R50X30R-C	ECG	Roll	50 x 30	Red	10
		R63-30R		R63-30R	ECG	Roll	63 x 30	Red	10
	SUZUKEN 302	R112X27Z		R112X27Z	ECG	Pack	112 x 90 x 300	Red	5
	SUZUKEN 301/106	R183X52Z-D		R183X52Z-D	ECG	Pack	183 x 270 x 200	Red	2
Kontron	STRESS TEST	204433		204433	ECG	Roll	216 x 45	Red	1
Macquarie	SPENGLER CT100	E5-003		E5-003	ECG	Roll	57 x 15	Green	5
Madison					AUD	Roll	60 x 10	Silver	100
Maico	MA 728/800/MA440				AUD	Roll	62 x 45	White	5
Marquette	MAC 6 / 12-15 / CASE 12-16	9402-020/034		9402-020/034	ECG	Pack	214 x 280 x 300	Red	8
		9402-046/026		9402-046/026	ECG	Roll	50 x 30	White	10
	MAC PC	9402-023		9402-023	ECG	Roll	108 x 23	Red	5
	MAC 12-15	9402-021		9402-021	ECG	Pack	210 x 280 x 300	Red	8
Mediana		YM4121		YM4121	ECG	Roll	108 x 25	Red	5
Medigate		Medigate		Medigate	ECG	Roll	111.5 x 25	Red	5
Mennen	1211-DG				ECG	Pack	220 x 300 x 500	White	1
	HORIZON 9000	914-159-010		914-159-010	ECG	Roll	50 x 30	Red	10
	CARDIOPAK 2000	800-040-060		800-040-060	ECG	Roll	63 x 30	Red	10
	HORIZON 1000/2000	800-040-010		800-040-010	ECG	Pack	216 x 280 x 500	White	1
M R L	CATHLAB	800-040-200		800-040-200	ECG	Roll	50 x 30	Red	10
		MR001729		MR001729	ECG	Roll	50 x 30	White	10
Micromedical	MICROLAB 3300	PSA 1600		PSA 1600	SPIRO	Roll	112 x 10	White	5
Mijnhart	VICATEST 2	226 084 02		226 084 02	SPIRO	Roll	60 x 25	White	5
	VICATEST 4/RESITEST	226 084 04		226 084 04	SPIRO	Roll	70 x 25	White	5
	VICATEST 3/5/P2	226 137 01		226 137 01	SPIRO	Roll	79 x 25	White	5
Minato	Microspiro	Hi 298/Hi498/AS500		Hi 298/Hi498/AS500	SPIRO	Roll	44 X 50OD	Bond	2
	Microspiro	Hi 198/AS303		Hi 198/AS303	SPIRO	Roll	55 x 50OD	Bond	5
	Autospiro	Hi 601/Hi501		Hi 601/Hi501	SPIRO	Roll	112 x 50OD	Bond	5
Mortara	ELI 50	00305003		00305003	ECG	Pack	50 x 70 x 200	Red	5
	ELI 100	00310801		00310801	ECG	Roll	108 x 23	Red	5
	ELI 150	9100-028-50		9100-028-50	ECG	Pack	108 x 140 x 214	Red	1
	ELI 200	00321504		00321504	ECG	Roll	215 x 50	Red	1
	ELI 200	00321505		00321505	ECG	Pack	215 x 280 x 200	Red	1
	ELI230	9100-029-50		9100-029-50	ECG	Roll	210 x 50OD	Red	1
	ELI 250	9100-026-51		9100-026-51	ECG	Pack	210 x 300 x 200	Red	1
	PORTAIT	00329701		00329701	ECG	Pack	280 x 215 x 500	Red	1
	ECG 1201 / 5201 / 5251 / 6511/6851	MS50-3		MS50-3	ECG	Roll	50 x 30	Green	10
	ECG 5303 / 5403	MS145-3		MS145-3	ECG	Roll	145 x 30	Orange	5
	ECG 6551	RQS50-3		RQS50-3	ECG	Roll	50 x 30	Orange	10
Nihon-Kohden	ECG 6651 / 6165	FQS50-3-75		FQS50-3-75	ECG	Pack	50 x 75 x 400	Orange	10
	TEC 7000SERIES/5500/5531/5521/7700/7721	FQS50-2-100		FQS50-2-100	ECG	Pack	50 x 100 x 200	Orange	10
	TEC 7521K/7100K/7511/7531/7621	FQS 50-3-100		FQS 50-3-100	ECG	Pack	50 x 100 x 300	Orange	10
	TEC7200/7100	FW 50-3-100		FW 50-3-100	ECG	Pack	50 x 100 x 300	White	10
	CARDIOFAX 8110 / ECG9620K	RQS 63-3		RQS 63-3	ECG	Roll	63 x 30	Orange	10

	Monitor / Recorder	Manufacturer #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Nihon-Kohden	ECG 8110	FQS 63-3-75	FQS 63-3-75	ECG	Pack	63 x 75 x 400	Orange	10
	ECG 6303 / 6403 / 6453	RQS145-3	RQS145-3	ECG	Roll	145 x 30	Orange	5
	ECG8400 / 8410 / 8420 / 8430	FQS 145-2.5-100	FQS 145-2.5-100	ECG	Pack	145 x 100 x 250	Orange	5
	ECG8240 / 8210 / 8220 / 8250	FQS 145-3A-150	FQS 145-3A-150	ECG	Pack	145 x 150 x 200	Orange	1
	ECG8370K	FQS216-3A-150	FQS216-3A-150	ECG	Pack	216 x 150 x 200	Orange	1
	ECG 8410 / 8420 / 8430 / 8210 / 8820	RQS 145-6B-150	RQS 145-6B-150	ECG	Roll	145 x 60	Orange	1
	ECG 8210 / 8220 / 8250 / 8270	RQS 145-6A-300	RQS 145-6A-300	ECG	Roll	145 x 60	Orange	1
	CTG-8100K	FQS150-2.7-90	FQS150-2.7-90	CTG	Pack	150 x 90 x 300	Green	5
	ECG8300 SERIES	RQS216-3B-150	RQS216-3B-150	ECG	Roll	216 x 30	Orange	1
	ECG 8300 / 8310 / 8320 / 8340	RQS216-6B-150	RQS216-6B-150	ECG	Roll	216 x 60	Orange	1
	ECG8300 / 8310 / 8320 / 8340	RQS216-6B-300	RQS216-6B-300	ECG	Roll	216 x 60	Orange	1
	ECG 6551	RQW50-3	RQW50-3	ECG	Roll	50 x 30	White	10
	9010K/9020K/1950/1250/9022	FQW 110-2-140	FQW 110-2-140	ECG	Pack	110 x 140 x 143	White	5
	NEUROPAK 2 / MEM 7102 K	FQW110-3-150	FQW110-3-150	ECG	Pack	110 x 150 x 200	White	5
	ECG9132/M1350/9130/9110/9230	FQW210-3-140	FQW210-3-140	ECG	Pack	210 x 140 x 215	White	1
		FQW210-3-295	FQW210-3-295	ECG	Pack	210 x 295 x 102	White	1
	MEM4104/WS180G/8320K	FQW 210-10A	FQW 210-10A	ECG	Pack	210 x 300 x 333	White	1
	ECG682G/800R	FQW210-30A	FQW 210-30A	ECG	Pack	210 x 300 x 1000	White	1
		FQW300-10A	FQW300-10A	EEG	Pack	300 x 300 x 333	White	1
		FQW210-7-295	FQW210-7-295	ECG	Pack	210 x 295 x 237	White	1
	ECG 9000/9320	FQW 210-10-295	FQW 210-10-295	ECG	Pack	210 x 300 x 333	White	1
	EEG7109.7209	FV 155 EEG	FV 155 EEG	EEG	Pack	155 x 300 x 1000	Green	1
	NEUROSCRIPT EE208	FV210 EEG	FV210 EEG	EEG	Pack	210 x 300 x 1000	Green	1
	TECA DG32	FV215 EEG	FV215 EEG	EEG	Pack	215 x 300 x 1000	Green	1
	EEG4414D.5214A.6714.7314	FV 245 EEG	FV 245 EEG	EEG	Pack	245 x 300 x 1000	Green	1
	EEG4414.4418.6714E.6718.4317.6514A	FV 345 EEG	FV 345 EEG	EEG	Pack	345 x 300 x 1000	Green	1
	EEG4221.4321.4421	FV 400 EEG	FV 400 EEG	EEG	Pack	400 x 300 x 1000	Green	1
		FV420 EEG	FV420 EEG	EEG	Pack	420 x 300 x 1000	Green	1
Nippon Collin	Refer to Collin							
Omron	Sphygs	864126	864126	BP	Roll	38 x 5	White	5
OTE Biomedical	REFER ESAOTE							
Oxford	Refer to Sonicaid							
Physio-Control	LIFE PAK 5	09-10419	09-10419	ECG	Roll	42 x 20	Red	5
	LIFE PAK 4/6/7/110	801-262	801-262	ECG	Roll	50 x 30	Red	10
	LIFE PAK 9/9P/10/300	804-700	804-700	ECG	Roll	50 x 30	Green	10
	LIFE PAK 11/112	805-319	805-319	ECG	Roll	108 x 25	Red	5
Philips	STRESS VUE	989803E +11	989803E +11	ECG	Pack	214 x 280 x 200	Red	1
Quinton	Q700	16631-002	16631-002	ECG	Pack	216 x 280 x 215	Orange	1
	Q750	18416-001	18416-001	ECG	Pack	216 x 280 x 215	Orange	1
	Q3000/Q4000	15022-001	15022-001	ECG	Pack	216 x 280 x 215	White	1
	Q4500	30146-001	30146-001	ECG	Pack	210 x 300 x 200	Orange	1
	Q5000	16211-002	16211-002	ECG	Pack	216 x 280 x 215	White	1

	Monitor / Recorder	Manufacturer #	Chart #	Re-order #	Category	Type Roll or Fanfold	Dimensions mm / m / OD / Sheets	Grid Colour	Packing Rolls / Packs Per Box
Siemens	CARDIOSTAT 701	43 55 053	43 55 053	43 55 053	ECG	Roll	60 x 45	Red	10
		42 55 053E3018	42 55 053E3018	42 55 053E3018	ECG	Roll	50 x 45	Red	10
	SIREDOC 220	45 26 208EH415	45 26 208EH415	45 26 208EH415	ECG	Pack	216 x 140 x 714	White	1
	CARDIOSTAT 31/31S	93 90 477E100E	93 90 477E100E	93 90 477E100E	ECG	Pack	104 x 100 x 300	Red	1
	CARDIOSTAT 31/31S	93 90 766E100E	93 90 766E100E	93 90 766E100E	ECG	Roll	104 x 30	Red	5
	MEGACART	61 74 903E100E	61 74 903E100E	61 74 903E100E	ECG	Pack	210 x 300 x 500	Orange	1
	SICARD 460	96 58 535E100E	96 58 535E100E	96 58 535E100E	ECG	Pack	147 x 100 x 400	Orange	2
Simonson & Weel	SICARD 440	96 20 741EH49E	96 20 741EH49E	96 20 741EH49E	ECG	Roll	147 x 25	Orange	5
	DMS 730 / 50	CP603	CP603	CP603	ECG	Roll	50 x 30	Yellow	10
	I Channel	CP1	CP1	CP1	ECG	Roll	63 x 45	Yellow	10
		CA999	CA999	CA999	ECG	Roll	50 x 25	White	5
	ATHENA	9541-01	9541-01	9541-01	ECG	Pack	120 x 150 x 300	White	2
Sonicaid	MERIDIAN / IMEX	8400/8003	8400/8003	8400/8003	CTG	Pack	143 x 150 x 300	White	2
Spacelabs	PC 90322	307568	307568	307568	ECG	Pack	120 x 150 x 300	White	2
		006-0196-00	006-0196-00	006-0196-00	ECG	Pack	50 x 120 x 150	White	1
	REPORT GENERATOR	90239A	90239A	90239A	ECG	Roll	57 x 10	White	5
Spirolife	201	VIS-C-004	VIS-C-004	VIS-C-0074	SPIRO	Roll	112 x 10	White	5
	220	MED-220	MED-220	MED-220	SPIRO	Roll	112 x 20	White	5
Spirometer		TR40	TR40	TR40	SPIRO	Roll	110 x 30	White	5
Spirometrics	FLOWMATE (Refer to Jaeger)								
Schwarzer		758201	758201	758201	ECG	Pack	215 x 280 x 215	White	1
Sunray		Sunray	Sunray		CTG	CTG	112x100x150	Orange	5
Technicon	Analyser	650-00-2601	650-00-2601	650-00-2601	LAB	Roll	57 x 30	White	4
	Analyser	650-00-2701	650-00-2701	650-00-2701	LAB	Roll	80 x 30	White	2
	Analyser	650-00-2801	650-00-2801	650-00-2801	LAB	Roll	110 x 30	White	2
Thymatron		ECF	ECF	ECF	ECG	Pack	108 x 150 x 300	White	1
Toitu	MT-430	0030-019	0030-019	0030-019	CTG	Pack	151 x 150 x 200	Green	5
		0030-022	0030-022	0030-022	CTG	Pack	151 x 150 x 200	Green	5
	MT516	0030-026	0030-026	0030-026	CTG	Pack	152 x 150 x 200	Green	1
	MT-325/MT-332MT320	0030-005	0030-005	0030-005	CTG	Pack	150 x 150 x 100	Green	5
Urillax	Urillax				LAB	Roll	58 x 30	White	5
Vitalograph	COMPAC	66149	66149	66149	SPIRO	Roll	112 x 25	White	5
		20.200	20.200	20.200	SPIRO	Sheet	198 x 218	Blue	200
	ALPHA	44058	44058	44058	SPIRO	Roll	112 x 35	White	5
	SPIROSIFT 500				SPIRO	Roll	110 x 25	White	5
Wakeling	FM210	9210-0006	9210-0006	9210-0006	ECG	Pack	100 x 100 x 150	Blue	5
	FM260	9260-0005	9260-0005	9260-0005	ECG	Pack	145 x 100 x 150	Blue	5
Weich Allyn	CP10-20	94016-0000	94016-0000	94016-0000	ECG	Pack	210 x 280 x 200	Red	1
	CP100/200	94018-0000	94018-0000	94018-0000	ECG	Pack	210 x 280 x 200	Red	1
	PEAK 20/30	1739	1739		ECG	Roll	50x30	White	10
Zoll	Zoll Series M	8000-0300	8000-0300	8000-0300	ECG	Pack	90 x 90 x 200	Red	5

Should you not find the product you are looking for listed in our catalogue please contact us.