


Phototherapy Unit

CHS-PU34




Functionality

- 3 steps brightness function
- Timer function
- Operating time and total usage time display



Adjustability

- Wide exposure range by tilting and height adjustment



Accessory

CONTENTS	EA
Eye Mask	Large 1, Small 1
Curtain	1
Manual	1
Option Accessory	EA
Lamp	1

Specification

CONTENTS	DESCRIPTION
Dimension(mm)	Lamp housing: 620(W) x 540(D) Base: 600(W) x 510(D) Height: 1,140~ 1,680
Lamp	Non ultra-violet ray lamp 4pcs
Power Source	AC 220V, 50/60Hz
Power Consumption	85W
Hour Meter	0 ~ 999hours
Caster	2"Swivel caster 4pcs
Housing Tilt	± 16°
Net Weight	25kg