

Digital LED Display Monitor type

HS-5020G / 4085G

Gravity Air Removal type _ 중력배기 방식



HS-5020G



HS-4085G

Specifications

MODEL		HS-5020G	HS-4085G
Style		Single Door type	
Chamber Volume		196 ℓ	85.4 ℓ
Operating Temperature		121.0℃ ~ 135.0℃	
Trouble Detection		Checked by Self Diagnosis Function	
Power Source	For Built-in Steam Source	AC 230V, 25A / AC 380V, 15A	AC 230V, 17A / AC 380V, 10A
	For External Steam Source	AC 230V, 5A	
Weight		350 kg	242 kg

HS-1636

Gravity Air Removal type _ 중력배기 방식

Military Portable Field Steam Sterilizer



Opened View

Closed View



HS-1636

Specifications

MODEL		HS-1636
Style		Portable for Field Use
Chamber Volume		118 ℓ (Cylindrical)
Operating Temperature		121.0℃ ~ 132.0℃
Power Source		AC 230V, 45A / AC 380V, 15A
Weight		155 kg