

# Large Steam Sterilizer

Touch Screen LCD Monitor type

## HS-5035PT / 5025PT / 5020PT / 4085PT

Automatic Sliding Door type\_ 수직 슬라이딩 자동도어



Single Door type



Double Door type



HS-5035PT



HS-5025PT



HS-5020PT



HS-4085PT

### Specifications

| MODEL                 |                           | HS-5035PT                            | HS-5025PT | HS-5020PT                   | HS-4085PT |
|-----------------------|---------------------------|--------------------------------------|-----------|-----------------------------|-----------|
| Style                 |                           | Automatic Vertical Sliding Door type |           |                             |           |
| Chamber Volume        |                           | 350 ℓ                                | 250 ℓ     | 196 ℓ                       | 85.4 ℓ    |
| Operating Temperature |                           | 121.0℃ ~ 135.0℃                      |           |                             |           |
| Trouble Detection     |                           | Checked by Self Diagnosis Function   |           |                             |           |
| Power Source          | For Built-in Steam Source | AC 230V, 55A / AC 380V, 32A          |           | AC 230V, 32A / AC 380V, 20A |           |
|                       | For External Steam Source | AC 230V, 5A                          |           |                             |           |
| Weight                | Single Door               | 803 kg                               | 710 kg    | 479 kg                      | 386 kg    |
|                       | Double Door               | 982 kg                               | 860 kg    | —                           | —         |