

Graphic LCD Monitor type

HS-5035 / 5025 / 5020 / 4085

Pre-vacuum & Gravity Air Removal type _ 선 진공 & 중력배기 방식



Single Door type



Double Door type



HS-5035



HS-5025



HS-5020



HS-4085

Specifications

MODEL		HS-5035	HS-5025	HS-5020	HS-4085
Style		Single Door type / Double Door type		Single Door type	
Chamber Volume		350 ℓ	250 ℓ	196 ℓ	85.4 ℓ
Operating Temperature		121.0℃ ~ 135.0℃			
Trouble Detection		Checked by Self Diagnosis Function			
Power Source	For Built-in Steam Source	AC 230V, 55A / AC 380V, 32A		AC 230V, 32A / AC 380V, 20A	AC 230V, 26A / AC 380V, 15A
	For External Steam Source	AC 230V, 5A			
Weight	Single Door	770 kg	665 kg	428 kg	334 kg
	Double Door	890 kg	755 kg		