

Graphic LCD Display Monitor type

HS-85SL / 60SL

First Vacuum Air Removal & Post-vacuum Drying type _ 진공배기 & 진공건조 방식



Automatic Door & Vertical Chamber type



HS-85SL



HS-60SL

Specifications

MODEL	HS-85SL	HS-60SL
Style	Vertical Chamber	
Chamber Volume	85 l (Cylindrical)	60 l (Cylindrical)
Operating Temperature	121.0℃ ~ 135.0℃	
Trouble Detection	Checked by Self Diagnosis Function	
Power Source	AC 230V / 3,300 W	
Weight	140 kg	130 kg

Digital LED Display Monitor type

HS-3460SD

First Vacuum Air Removal & Post-vacuum Drying type _ 진공배기 & 진공건조 방식



HS-3460SD

Specifications

MODEL	HS-3460SD
Style	Horizontal Chamber
Chamber Volume	60 l (Cylindrical)
Operating Temperature	121.0℃ ~ 135.0℃
Trouble Detection	Checked by Self Diagnosis Function
Power Source	AC 230V / 2,950 W
Weight	152 kg